




FRIZKON
INTERNATIONAL



LIPOSUCTION INSTRUMENTS

www.frizkontools.com
info@frizkontools.com



FRIZKON TOOLS

Established in 2006 FRIZKON TOOLS provides premium solutions for exceptional value to individuals and businesses. We are the third generation Ophthalmic, Electro Surgical and Plastic Surgery Instruments manufacturing company. Frizkon portfolio encompasses an experience of over 15 years to introduce technological advancements through wide range of product, we gather our strength from creative thought processing, strategic marketing techniques and experience of a team of individuals talented by nature and trained by the profession.

Position ourselves to be regional players of quality making products, improve on contemporary measure including cost, quality, service, speed of delivery and mobilization. We are customer driven we go the extra mile to make sure our clients exceptions are met and exceeded on every issue. Experience has shown that most clients prefer to work with our company rather than having many suppliers / companies. To achieve this, Frizkon have brought all your technological needs under one roof, we make sure that all our processes and methods conform to the highest ideals of professional behavior.

The commitment to excellence make us a world class well known company for our high quality standards from design to finished products and from sales to shipment all procedures and processes are organized to give our customers reliable products and a convenient and punctual services.

After years of expertise and experience, we still remain dedicated to continuous business and personal development, building strong customer relations and investing in research and development, Innovative ideas and personal dedication combined with high quality products is our specialty, our commitment of introducing quality products combined with honest, reliable service has proven to be a successful theory in ensuring the continued rapid growth of the company. Our transparent business policies together with highly professional and values are sheet anchor of our company management.

FEATURES OF OUR QUALITY CONTROL PROCEDURES

1:- MATERIALS

Choosing material is the main factor in the manufacturing process normally best kind of stain less is used in manufacturing the items, but other kind of French and Japanese steel and a series of 300 / 410 and 420 are also used as per the customer requirements PQA is the competent authority for the testing of raw materials, A rigid quality control policy is being adopted before releasing to the production process. Material test is to check the standard of steel as per desired specification and requirement of the customers.

2:- FORGING AND HAMMERING

Hot and cold forging is performed to the raw material (stainless steel bars) for the standard press operations in forming deep drawing of all type of instruments thru the specialized craftsmen. They annealed the forging process and prepare the materials for milling process.

3:- MILLING PROCESS.

Milling process is performed by the highly skilled machinist they are fully trained in accurately forming the male-female portions of box-joints forceps as well as serrations- retches and other wide range of instruments in uniformity. Almost all types of mechanical works are performed at the same place.

4:- FILLING.

Hand-made fabrication makes almost all items precise and finest in look and that is done by specialized and skilled persons who assemble male /female portions after adjusting of extracting forceps and scissors, tips, serrations, retches and all other type of items in uniformity and file precisely to make it semi-finish under strict quality control specifications.

5:-TEMPERING

All instruments are being tempered maintaining "Rock-Well Hardness". The tempering varies from 45 to 55 for various types of instruments. A hardness test is taken after tempering of each lot to maintain the uniformity. Further Passivation or a chemical process is performed to clean the items chemically.

6:- RIVETING

A wide range of assembly operation is performed from prototype quantities to production run. Each instrument is riveted / welded or spot welded precisely according to the requirement.

7:- GRINDING / POLISHING

Finally, to get the best finished products it is necessary to facilitate the craftsman with précised grinding and polishing equipment's for high uniformity utilizing best of polishing luster. Electro-plating and glass head process is used for the shined and matt finished steel instruments.

We are proud to say that Let us serve you to the best of our abilities.

For any query, call or mail us: info@frizkontools.com or visit www.frizkontools.com

FRIZKON INTERNATIONAL

Contact:- +92 52 3306091

+92 3325334864

+92 3436467358

Table of Contents

- A. Liposuction Cannulas
 - a) Infiltration Cannulas
 - b) General Suction Cannulas
 - c) Harvesting Cannulas
 - e) Injecting Cannulas
 - f) Speciality Cannulas

- B. Accessories

- C. Our Specialized Sets

Cannula Hubs/ Connections



Luer lock



Luer lock
(Stainless Steel)



Ultra Luer lock



Treaded



60cc Tommey hub

Machine Connections



FT-1000.01

Cannulas and Harvesters as a complete instrument ready to use with an aspiration machine for liposuction and large volume fat harvesting procedures.

Fixed Handle

Luer Lock Hub/Connector:

The universal Luer Lock hub enables cannula to attach directly with any Luer Lock Syringe.

Ultra Luer Lock Hub/Connector:

A modification to general Luer Lock hub which is more powerful, more stable, Leak proof, free from instability and contamination.

Threaded Hub/Connector:

Threaded connector hub style is used mostly for machine assisted Liposuction cannulas, FRIZKON INTERNATIONAL unique Light weight threaded cannula can also be used with ordinary syringe. All you need is FRIZKON INTERNATIONAL connector from threaded to Luer style.

60cc Tommey Hub/Connector:

An ideal hub style for surgeons why love to handle large fat volumes. This hub style connect cannula to a 60cc Tommey to syringe.


Fixed Handle Cannulas:

FRIZKON INTERNATIONAL super power grip handle fixed with the cannula of your own choice! The unique design of handle offer powerful grip harvesting large volumes. You will feel amazingly comfortable ever after a long operational time.

The thumb grip provide a unique opportunity to control the thrust.

INFILTRATION CANNULAS

Frizkon International Infiltration Cannulas are designed for safe and efficient delivery of tumescent solution. The spiral micro holes on cannula ensure an even distribution of local anesthetic which results in a longer anesthetic effect, less pain and a more pleasant liposuction experience for the patient.



Tumescent Infiltrator

This cannula is used specifically for the iRlusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60cc Tommeyhub

FT-101.01	Ø 2.5 mm X 15 cm
FT-101.02	Ø 2.5 mm X 20 cm
FT-101.03	Ø 3.0 mm X 20 cm
FT-101.04	Ø 3.0 mm X 25 cm

Treaded Hub

FT-102.01	Ø 2.0 mm X 15 cm
FT-102.02	Ø 3.0 mm X 15 cm
FT-102.03	Ø 3.0 mm X 20 cm
FT-102.04	Ø 3.0 mm X 30 cm
FT-102.05	Ø 4.0 mm X 20 cm
FT-102.06	Ø 4.0 mm X 30 cm

Ultra Luer lock

FT-103.01	Ø 1.0 mm X 5 cm
FT-103.02	Ø 1.0 mm X 7 cm
FT-103.03	Ø 1.5 mm X 7 cm
FT-103.04	Ø 1.5 mm X 9 cm

Luer lock

FT-104.01	Ø 1.0 mm X 4 cm
FT-104.02	Ø 1.0 mm X 5 cm
FT-104.03	Ø 1.5 mm X 7 cm
FT-104.04	Ø 1.5 mm X 9 cm
FT-104.05	Ø 1.5 mm X 10 cm
FT-104.06	Ø 1.5 mm X 12 cm
FT-104.07	Ø 2.0 mm X 10 cm
FT-104.08	Ø 2.0 mm X 12 cm
FT-104.09	Ø 2.0 mm X 15 cm
FT-104.10	Ø 2.0 mm X 20 cm
FT-104.11	Ø 2.5 mm X 15 cm
FT-104.12	Ø 2.5 mm X 20 cm
FT-104.13	Ø 3.0 mm X 15 cm
FT-104.14	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

FT-105.01	Ø 1.0 mm X 4 cm
FT-105.02	Ø 1.0 mm X 5 cm
FT-105.03	Ø 1.0 mm X 7 cm
FT-105.04	Ø 1.5 mm X 7 cm
FT-105.05	Ø 1.5 mm X 9 cm
FT-105.06	Ø 1.5 mm X 10 cm
FT-105.07	Ø 1.5 mm X 12 cm
FT-105.08	Ø 2.0 mm X 10 cm
FT-105.09	Ø 2.0 mm X 12 cm
FT-105.10	Ø 2.0 mm X 15 cm
FT-105.11	Ø 2.0 mm X 20 cm
FT-105.12	Ø 2.5 mm X 15 cm
FT-105.13	Ø 2.5 mm X 20 cm
FT-105.14	Ø 3.0 mm X 15 cm
FT-105.15	Ø 3.0 mm X 20 cm

Power Cannula

FT-106.01	Ø 2.0 mm X 15 cm
FT-106.02	Ø 2.0 mm X 20 cm
FT-106.03	Ø 2.0 mm X 25 cm
FT-106.04	Ø 2.0 mm X 30 cm
FT-106.05	Ø 2.0 mm X 35 cm
FT-106.06	Ø 3.0 mm X 15 cm
FT-106.07	Ø 3.0 mm X 20 cm
FT-106.08	Ø 3.0 mm X 25 cm
FT-106.09	Ø 3.0 mm X 30 cm
FT-106.10	Ø 3.0 mm X 35 cm
FT-106.11	Ø 4.0 mm X 15 cm
FT-106.12	Ø 4.0 mm X 20 cm
FT-106.13	Ø 4.0 mm X 25 cm
FT-106.14	Ø 4.0 mm X 30 cm
FT-106.15	Ø 4.0 mm X 35 cm
FT-106.16	Ø 5.0 mm X 15 cm
FT-106.17	Ø 5.0 mm X 20 cm
FT-106.18	Ø 5.0 mm X 25 cm
FT-106.19	Ø 5.0 mm X 30 cm
FT-106.20	Ø 5.0 mm X 35 cm
FT-106.21	Ø 6.0 mm X 15 cm
FT-106.22	Ø 6.0 mm X 20 cm
FT-106.23	Ø 6.0 mm X 25 cm
FT-106.24	Ø 6.0 mm X 30 cm



Tumescent Bayonet Infiltrator

The bayonet shaped handles allow an easier manipulation

60cc Tommeyhub

FT-101.11	Ø 2.0 mm X 15 cm
FT-101.12	Ø 2.5 mm X 20 cm
FT-101.13	Ø 3.0 mm X 20 cm
FT-101.14	Ø 3.0 mm X 25 cm

Treaded Hub

FT-102.11	Ø 2.0mm X 15 cm
FT-102.12	Ø 3.0mm X 15 cm
FT-102.13	Ø 3.0mm X 20 cm
FT-102.14	Ø 3.0mm X 30 cm
FT-102.15	Ø 4.0mm X 20 cm
FT-102.16	Ø 4.0mm X 30 cm

Ultra Luer lock

FT-103.11	Ø 1.0 mm X 5 cm
FT-103.12	Ø 1.0 mm X 7 cm
FT-103.13	Ø 1.5 mm X 7 cm
FT-103.14	Ø 1.5 mm X 9 cm

Luer lock


FT-104.11	Ø 1.0 mm X 4 cm
FT-104.12	Ø 1.0 mm X 5 cm
FT-104.13	Ø 1.0 mm X 7 cm
FT-104.14	Ø 1.5 mm X 7 cm
FT-104.15	Ø 1.5 mm X 9 cm
FT-104.16	Ø 1.5 mm X 10 cm
FT-104.17	Ø 1.5 mm X 12 cm
FT-104.18	Ø 2.0 mm X 10 cm
FT-104.19	Ø 2.0 mm X 12 cm
FT-104.20	Ø 2.0 mm X 15 cm
FT-104.21	Ø 2.0 mm X 20 cm
FT-104.22	Ø 2.5 mm X 15 cm
FT-104.23	Ø 2.5 mm X 20 cm
FT-104.24	Ø 3.0 mm X 15 cm
FT-104.25	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

FT-105.21	Ø 1.0 mm X 4 cm
FT-105.22	Ø 1.0 mm X 5 cm
FT-105.23	Ø 1.0 mm X 7 cm
FT-105.24	Ø 1.5 mm X 7 cm
FT-105.25	Ø 1.5 mm X 9 cm
FT-105.26	Ø 1.5 mm X 10 cm
FT-105.27	Ø 1.5 mm X 12 cm
FT-105.28	Ø 2.0 mm X 10 cm
FT-105.29	Ø 2.0 mm X 12 cm
FT-105.30	Ø 2.0 mm X 15 cm
FT-105.31	Ø 2.0 mm X 20 cm
FT-105.32	Ø 2.5 mm X 15 cm
FT-105.33	Ø 2.5 mm X 20 cm
FT-105.34	Ø 3.0 mm X 15 cm
FT-105.35	Ø 3.0 mm X 20 cm

Power Cannula

FT-106.31	Ø 2 mm X 15 cm
FT-106.32	Ø 2 mm X 20 cm
FT-106.33	Ø 2 mm X 25 cm
FT-106.34	Ø 2 mm X 30 cm
FT-106.35	Ø 2 mm X 35 cm
FT-106.36	Ø 3 mm X 15 cm
FT-106.37	Ø 3 mm X 20 cm
FT-106.38	Ø 3 mm X 25 cm
FT-106.39	Ø 3 mm X 30 cm
FT-106.40	Ø 3 mm X 35 cm
FT-106.41	Ø 4 mm X 15 cm
FT-106.42	Ø 4 mm X 20 cm
FT-106.43	Ø 4 mm X 25 cm
FT-106.44	Ø 4 mm X 30 cm
FT-106.45	Ø 4 mm X 35 cm
FT-106.46	Ø 5 mm X 15 cm
FT-106.47	Ø 5 mm X 20 cm
FT-106.48	Ø 5 mm X 25 cm
FT-106.49	Ø 5 mm X 30 cm
FT-106.50	Ø 5 mm X 35 cm
FT-106.51	Ø 6 mm X 15 cm
FT-106.52	Ø 6 mm X 20 cm
FT-106.53	Ø 6 mm X 25 cm
FT-106.54	Ø 6 mm X 30 cm
FT-106.55	Ø 6 mm X 35 cm



Tumescent Infiltrator Angled

This cannula is used specifically for the iRLusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60cc Tommeyhub

FT-101.21	Ø 2.0 mm X 15 cm
FT-101.22	Ø 2.5 mm X 20 cm
FT-101.23	Ø 3.0 mm X 20 cm
FT-101.24	Ø 3.0 mm X 25 cm

Treaded Hub

FT-102.21	Ø 2.0mm X 15 cm
FT-102.22	Ø 3.0mm X 15 cm
FT-102.23	Ø 3.0mm X 20 cm
FT-102.24	Ø 3.0mm X 30 cm
FT-102.25	Ø 4.0mm X 20 cm
FT-102.26	Ø 4.0mm X 30 cm

Ultra Luer lock

FT-103.21	Ø 1.0 mm X 5 cm
FT-103.22	Ø 1.0 mm X 7 cm
FT-103.23	Ø 1.5 mm X 7 cm
FT-103.24	Ø 1.5 mm X 9 cm

Luer lock


FT-104.31	Ø 1.0 mm X 4 cm
FT-104.32	Ø 1.0 mm X 5 cm
FT-104.33	Ø 1.0 mm X 7 cm
FT-104.34	Ø 1.5 mm X 7 cm
FT-104.35	Ø 1.5 mm X 9 cm
FT-104.36	Ø 1.5 mm X 10 cm
FT-104.37	Ø 1.5 mm X 12 cm
FT-104.38	Ø 2.0 mm X 10 cm
FT-104.39	Ø 2.0 mm X 12 cm
FT-104.40	Ø 2.0 mm X 15 cm
FT-104.41	Ø 2.0 mm X 20 cm
FT-104.42	Ø 2.5 mm X 15 cm
FT-104.43	Ø 2.5 mm X 20 cm
FT-104.44	Ø 3.0 mm X 15 cm
FT-104.45	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

FT-105.41	Ø 1.0 mm X 4 cm
FT-105.42	Ø 1.0 mm X 5 cm
FT-105.43	Ø 1.0 mm X 7 cm
FT-105.44	Ø 1.5 mm X 7 cm
FT-105.45	Ø 1.5 mm X 9 cm
FT-105.46	Ø 1.5 mm X 10 cm
FT-105.47	Ø 1.5 mm X 12 cm
FT-105.48	Ø 2.0 mm X 10 cm
FT-105.49	Ø 2.0 mm X 12 cm
FT-105.50	Ø 2.0 mm X 15 cm
FT-105.51	Ø 2.0 mm X 20 cm
FT-105.52	Ø 2.5 mm X 15 cm
FT-105.53	Ø 2.5 mm X 20 cm
FT-105.54	Ø 3.0 mm X 15 cm
FT-105.55	Ø 3.0 mm X 20 cm

Power Cannula

FT-106.61	Ø 2 mm X 15 cm
FT-106.62	Ø 2 mm X 20 cm
FT-106.63	Ø 2 mm X 25 cm
FT-106.64	Ø 2 mm X 30 cm
FT-106.65	Ø 2 mm X 35 cm
FT-106.66	Ø 3 mm X 15 cm
FT-106.67	Ø 3 mm X 20 cm
FT-106.68	Ø 3 mm X 25 cm
FT-106.69	Ø 3 mm X 30 cm
FT-106.70	Ø 3 mm X 35 cm
FT-106.71	Ø 4 mm X 15 cm
FT-106.72	Ø 4 mm X 20 cm
FT-106.73	Ø 4 mm X 25 cm
FT-106.74	Ø 4 mm X 30 cm
FT-106.75	Ø 4 mm X 35 cm
FT-106.76	Ø 5 mm X 15 cm
FT-106.77	Ø 5 mm X 20 cm
FT-106.78	Ø 5 mm X 25 cm
FT-106.79	Ø 5 mm X 30 cm
FT-106.80	Ø 5 mm X 35 cm
FT-106.81	Ø 6 mm X 15 cm
FT-106.82	Ø 6 mm X 20 cm
FT-106.83	Ø 6 mm X 25 cm
FT-106.84	Ø 6 mm X 30 cm



Facelift Infiltration Cannula

This cannula is used specifically for the iRlusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60ccTommeyhub

FT-101.31	Ø 2.0 mm X 15 cm
FT-101.32	Ø 2.5 mm X 20 cm
FT-101.33	Ø 3.0 mm X 20 cm
FT-101.34	Ø 3.0 mm X 25 cm

Treaded Hub

FT-102.31	Ø 2.0mm X 15 cm
FT-102.32	Ø 3.0mm X 15 cm
FT-102.33	Ø 3.0mm X 20 cm
FT-102.34	Ø 3.0mm X 30 cm
FT-102.35	Ø 4.0mm X 20 cm
FT-102.36	Ø 4.0mm X 30 cm

Ultra Luer lock

FT-103.31	Ø 1.0 mm X 5 cm
FT-103.32	Ø 1.0 mm X 7 cm
FT-103.33	Ø 1.5 mm X 7 cm
FT-103.34	Ø 1.5 mm X 9 cm

Luer lock

FT-104.51	Ø 1.0 mm X 4 cm
FT-104.52	Ø 1.0 mm X 5 cm
FT-104.53	Ø 1.0 mm X 7 cm
FT-104.54	Ø 1.5 mm X 7 cm
FT-104.55	Ø 1.5 mm X 9 cm
FT-104.56	Ø 1.5 mm X 10 cm
FT-104.57	Ø 1.5 mm X 12 cm
FT-104.58	Ø 2.0 mm X 10 cm
FT-104.59	Ø 2.0 mm X 12 cm
FT-104.50	Ø 2.0 mm X 15 cm
FT-104.51	Ø 2.0 mm X 20 cm
FT-104.52	Ø 2.5 mm X 15 cm
FT-104.53	Ø 2.5 mm X 20 cm
FT-104.54	Ø 3.0 mm X 15 cm
FT-104.55	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

FT-105.61	Ø 1.0 mm X 4 cm
FT-105.62	Ø 1.0 mm X 5 cm
FT-105.63	Ø 1.0 mm X 7 cm
FT-105.64	Ø 1.5 mm X 7 cm
FT-105.65	Ø 1.5 mm X 9 cm
FT-105.66	Ø 1.5 mm X 10 cm
FT-105.67	Ø 1.5 mm X 12 cm
FT-105.68	Ø 2.0 mm X 10 cm
FT-105.69	Ø 2.0 mm X 12 cm
FT-105.70	Ø 2.0 mm X 15 cm
FT-105.71	Ø 2.0 mm X 20 cm
FT-105.72	Ø 2.5 mm X 15 cm
FT-105.73	Ø 2.5 mm X 20 cm
FT-105.74	Ø 3.0 mm X 15 cm
FT-105.75	Ø 3.0 mm X 20 cm

Power Cannula

FT-106.91	Ø 2 mm X 15 cm
FT-106.92	Ø 2 mm X 20 cm
FT-106.93	Ø 2 mm X 25 cm
FT-106.94	Ø 2 mm X 30 cm
FT-106.95	Ø 2 mm X 35 cm
FT-106.96	Ø 3 mm X 15 cm
FT-106.97	Ø 3 mm X 20 cm
FT-106.98	Ø 3 mm X 25 cm
FT-106.99	Ø 3 mm X 30 cm
FT-106.100	Ø 3 mm X 35 cm
FT-106.101	Ø 4 mm X 15 cm
FT-106.102	Ø 4 mm X 20 cm
FT-106.103	Ø 4 mm X 25 cm
FT-106.104	Ø 4 mm X 30 cm
FT-106.105	Ø 4 mm X 35 cm
FT-106.106	Ø 5 mm X 15 cm
FT-106.107	Ø 5 mm X 20 cm
FT-106.108	Ø 5 mm X 25 cm
FT-106.109	Ø 5 mm X 30 cm
FT-106.110	Ø 5 mm X 35 cm
FT-106.111	Ø 6 mm X 15 cm
FT-106.112	Ø 6 mm X 20 cm
FT-106.113	Ø 6 mm X 25 cm
FT-106.114	Ø 6 mm X 30 cm
FT-106.115	Ø 6 mm X 35 cm



GENERAL SUCTION CANNULA

Frizkon International General suction cannula are ideal tool for large volume Harvesting with ease perfection.

These cannulas are designed to facilitate you to get perfectly designed goals in minimum time.



StevensStyleCannula

60ccTommeyhub

FT-201.01	Ø 2.0 mm X 15 cm
FT-201.02	Ø 2.5 mm X 20 cm
FT-201.03	Ø 3.0 mm X 20 cm
FT-201.04	Ø 3.0 mm X 25 cm

Treaded Hub

FT-202.01	Ø 3.0 mm X 20 cm
FT-202.02	Ø 3.0 mm X 30 cm
FT-202.03	Ø 4.0 mm X 30 cm
FT-202.04	Ø 5.0 mm X 30 cm
FT-202.05	Ø 6.0 mm X 30 cm

Ultra Luer lock

FT-203.01	Ø 1.0 mm X 5 cm
FT-203.02	Ø 1.0 mm X 7 cm
FT-203.03	Ø 1.5 mm X 7 cm
FT-203.04	Ø 1.5 mm X 9 cm
FT-203.05	Ø 1.5 mm X 7 cm

Luer lock


FT-204.01	Ø 1.0 mm X 4 cm
FT-204.02	Ø 1.0 mm X 5 cm
FT-204.03	Ø 1.0 mm X 7 cm
FT-204.04	Ø 1.5 mm X 7 cm
FT-204.05	Ø 1.5 mm X 9 cm
FT-204.06	Ø 1.5 mm X 10 cm
FT-204.07	Ø 1.5 mm X 12 cm
FT-204.08	Ø 2.0 mm X 10 cm
FT-204.09	Ø 2.0 mm X 12 cm
FT-204.10	Ø 2.0 mm X 15 cm
FT-204.11	Ø 2.0 mm X 20 cm
FT-204.12	Ø 2.5 mm X 15 cm
FT-204.13	Ø 2.5 mm X 20 cm
FT-204.14	Ø 3.0 mm X 15 cm
FT-204.15	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

FT-205.01	Ø 1.0 mm X 4 cm
FT-205.02	Ø 1.0 mm X 5 cm
FT-205.03	Ø 1.0 mm X 7 cm
FT-205.04	Ø 1.5 mm X 7 cm
FT-205.05	Ø 1.5 mm X 9 cm
FT-205.06	Ø 1.5 mm X 10 cm
FT-205.07	Ø 1.5 mm X 12 cm
FT-205.08	Ø 2.0 mm X 10 cm
FT-205.09	Ø 2.0 mm X 12 cm
FT-205.10	Ø 2.0 mm X 15 cm
FT-205.11	Ø 2.0 mm X 20 cm
FT-205.12	Ø 2.5 mm X 15 cm
FT-205.13	Ø 2.5 mm X 20 cm
FT-205.14	Ø 3.0 mm X 15 cm
FT-205.15	Ø 3.0 mm X 20 cm

Power Cannula

FT-206.01	Ø 2 mm X 15 cm
FT-206.02	Ø 2 mm X 20 cm
FT-206.03	Ø 2 mm X 25 cm
FT-206.04	Ø 2 mm X 30 cm
FT-206.05	Ø 2 mm X 35 cm
FT-206.06	Ø 3 mm X 15 cm
FT-206.07	Ø 3 mm X 20 cm
FT-206.08	Ø 3 mm X 25 cm
FT-206.09	Ø 3 mm X 30 cm
FT-206.10	Ø 3 mm X 35 cm
FT-206.11	Ø 4 mm X 15 cm
FT-206.12	Ø 4 mm X 20 cm
FT-206.13	Ø 4 mm X 25 cm
FT-206.14	Ø 4 mm X 30 cm
FT-206.15	Ø 4 mm X 35 cm
FT-206.16	Ø 5 mm X 15 cm
FT-206.17	Ø 5 mm X 20 cm
FT-206.18	Ø 5 mm X 25 cm
FT-206.19	Ø 5 mm X 30 cm
FT-206.20	Ø 5 mm X 35 cm
FT-206.21	Ø 6 mm X 15 cm
FT-206.22	Ø 6 mm X 20 cm
FT-206.23	Ø 6 mm X 25 cm
FT-206.24	Ø 6 mm X 30 cm
FT-206.25	Ø 6 mm X 35 cm



Cobra Bibevel Style

More aggressive than the Standard, Pyramid or Mercedes, this cannula has a tapered tip with two ports and works effectively within dense fibrous tissue.

60cc Tommeyhub

FT-201.11	Ø 2.5 mm X 15 cm5
FT-201.12	Ø 2.5 mm X 20 cm
FT-201.13	Ø 3.0 mm X 15 cm
FT-201.14	Ø 3.0 mm X 20 cm
FT-201.15	Ø 3.0 mm X 25 cm
FT-201.16	Ø 4.0 mm X 25 cm
FT-201.17	Ø 4.0 mm X 30 cm
FT-201.18	Ø 4.0 mm X 35 cm
FT-201.19	Ø 4.5 mm X 30 cm
FT-201.20	Ø 4.5 mm X 35 cm

Treaded Hub

FT-202.11	Ø 3.0 mm X 15 cm
FT-202.12	Ø 3.0 mm X 25 cm
FT-202.13	Ø 4.0 mm X 15 cm
FT-202.14	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.11	Ø 2.0 mm X 15 cm
FT-203.12	Ø 2.5 mm X 15 cm
FT-203.13	Ø 2.5 mm X 20 cm
FT-203.14	Ø 3.0 mm X 15 cm
FT-203.15	Ø 3.0 mm X 20 cm
FT-203.16	Ø 3.0 mm X 25 cm

Luer lock


FT-204.21	Ø 2.0 mm X 15 cm
FT-204.22	Ø 2.5 mm X 15 cm
FT-204.23	Ø 2.5 mm X 20 cm
FT-204.24	Ø 3.0 mm X 15 cm
FT-204.25	Ø 3.0 mm X 20 cm
FT-204.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-205.21	Ø 2.0 mm X 15 cm
FT-205.22	Ø 2.5 mm X 15 cm
FT-205.23	Ø 2.5 mm X 20 cm
FT-205.24	Ø 3.0 mm X 15 cm
FT-205.25	Ø 3.0 mm X 20 cm
FT-205.26	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.31	Ø 2mm X 15 cm
FT-206.32	Ø 2mm X 20 cm
FT-206.33	Ø 2mm X 25 cm
FT-206.34	Ø 2mm X 30 cm
FT-206.35	Ø 2mm X 35 cm
FT-206.36	Ø 3mm X 15 cm
FT-206.37	Ø 3mm X 20 cm
FT-206.38	Ø 3mm X 25 cm
FT-206.39	Ø 3mm X 30 cm
FT-206.40	Ø 3mm X 35 cm
FT-206.41	Ø 4mm X 15 cm
FT-206.42	Ø 4mm X 20 cm
FT-206.43	Ø 4mm X 25 cm
FT-206.44	Ø 4mm X 30 cm
FT-206.45	Ø 4mm X 35 cm
FT-206.46	Ø 5mm X 15 cm
FT-206.47	Ø 5mm X 20 cm
FT-206.48	Ø 5mm X 25 cm
FT-206.49	Ø 5mm X 30 cm
FT-206.50	Ø 5mm X 35 cm
FT-206.51	Ø 6mm X 15 cm
FT-206.52	Ø 6mm X 20 cm
FT-206.53	Ø 6mm X 25 cm
FT-206.54	Ø 6mm X 30 cm
FT-206.55	Ø 6mm X 35 cm



Cobra Style Round Tip

Less aggressive than the Cobra Bibevel this cell friendly cannula has a rounded tip with two ports and works effectively within dense fibrous tissue during liposuction procedures.

60ccTommeyhub

FT-201.21	Ø 2.5 mm X 15 cm
FT-201.22	Ø 2.5 mm X 20 cm
FT-201.23	Ø 3.0 mm X 15 cm
FT-201.24	Ø 3.0 mm X 20 cm
FT-201.25	Ø 3.0 mm X 25 cm
FT-201.26	Ø 4.0 mm X 25 cm
FT-201.27	Ø 4.0 mm X 30 cm
FT-201.28	Ø 4.0 mm X 35 cm
FT-201.29	Ø 4.5 mm X 30 cm
FT-201.30	Ø 4.5 mm X 35 cm

Treaded Hub

FT-202.21	Ø 3.0 mm X 15 cm
FT-202.22	Ø 3.0 mm X 25 cm
FT-202.23	Ø 4.0 mm X 15 cm
FT-202.24	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.21	Ø 2.0 mm X 15 cm
FT-203.22	Ø 2.5 mm X 15 cm
FT-203.23	Ø 2.5 mm X 20 cm
FT-203.24	Ø 3.0 mm X 15 cm
FT-203.25	Ø 3.0 mm X 20 cm
FT-203.26	Ø 3.0 mm X 25 cm

Luer lock

FT-204.31	Ø 2.0 mm X 15 cm
FT-204.32	Ø 2.5 mm X 15 cm
FT-204.33	Ø 2.5 mm X 20 cm
FT-204.34	Ø 3.0 mm X 15 cm
FT-204.35	Ø 3.0 mm X 20 cm
FT-204.36	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-205.31	Ø 2.0 mm X 15 cm
FT-205.32	Ø 2.5 mm X 15 cm
FT-205.33	Ø 2.5 mm X 20 cm
FT-205.34	Ø 3.0 mm X 15 cm
FT-205.35	Ø 3.0 mm X 20 cm
FT-205.36	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.61	Ø 2 mm X 15 cm
FT-206.62	Ø 2 mm X 20 cm
FT-206.63	Ø 2 mm X 25 cm
FT-206.64	Ø 2 mm X 30 cm
FT-206.65	Ø 2 mm X 35 cm
FT-206.66	Ø 3 mm X 15 cm
FT-206.67	Ø 3 mm X 20 cm
FT-206.68	Ø 3 mm X 25 cm
FT-206.69	Ø 3 mm X 30 cm
FT-206.70	Ø 3 mm X 35 cm
FT-206.71	Ø 4 mm X 15 cm
FT-206.72	Ø 4 mm X 20 cm
FT-206.73	Ø 4 mm X 25 cm
FT-206.74	Ø 4 mm X 30 cm
FT-206.75	Ø 4 mm X 35 cm
FT-206.76	Ø 5 mm X 15 cm
FT-206.77	Ø 5 mm X 20 cm
FT-206.78	Ø 5 mm X 25 cm
FT-206.79	Ø 5 mm X 30 cm
FT-206.80	Ø 5 mm X 35 cm
FT-206.81	Ø 6 mm X 15 cm
FT-206.82	Ø 6 mm X 20 cm
FT-206.83	Ø 6 mm X 25 cm
FT-206.84	Ø 6 mm X 30 cm
FT-206.85	Ø 6 mm X 35 cm



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula



Standard, Single Hole

This is the most common cell friendly cannula design for soft tissue removal. The single port is located 180 degrees from the thumb rest on the hub.

60cc Tommeyhub

FT-201.31	Ø 2.5 mm X 15 cm
FT-201.32	Ø 2.5 mm X 20 cm
FT-201.33	Ø 3.0 mm X 15 cm
FT-201.34	Ø 3.0 mm X 20 cm
FT-201.35	Ø 3.0 mm X 25 cm
FT-201.36	Ø 4.0 mm X 25 cm
FT-201.37	Ø 4.0 mm X 30 cm
FT-201.38	Ø 4.0 mm X 35 cm
FT-201.39	Ø 4.5 mm X 30 cm
FT-201.40	Ø 4.5 mm X 35 cm

Luer lock (Stainless Steel)

FT-205.41	Ø 1.5 mm X 12 cm
FT-205.42	Ø 2.0 mm X 15 cm
FT-205.43	Ø 2.0 mm X 20 cm
FT-205.44	Ø 2.5 mm X 15 cm
FT-205.45	Ø 2.5 mm X 20 cm
FT-205.46	Ø 3.0 mm X 15 cm
FT-205.47	Ø 3.0 mm X 20 cm
FT-205.48	Ø 3.0 mm X 25 cm

Treaded Hub

FT-202.31	Ø 3.0 mm X 15 cm
FT-202.32	Ø 3.0 mm X 25 cm
FT-202.33	Ø 4.0 mm X 15 cm
FT-202.34	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.31	Ø 1.5 mm X 12 cm
FT-203.32	Ø 2.0 mm X 15 cm
FT-203.33	Ø 2.0 mm X 20 cm
FT-203.34	Ø 2.5 mm X 15 cm
FT-203.35	Ø 2.5 mm X 20 cm
FT-203.36	Ø 3.0 mm X 15 cm
FT-203.37	Ø 3.0 mm X 20 cm
FT-203.38	Ø 3.0 mm X 25 cm

Luer lock

FT-204.41	Ø 1.5 mm X 12 cm
FT-204.42	Ø 2.0 mm X 15 cm
FT-204.43	Ø 2.0 mm X 20 cm
FT-204.44	Ø 2.5 mm X 15 cm
FT-204.45	Ø 2.5 mm X 20 cm
FT-204.46	Ø 3.0 mm X 15 cm
FT-204.47	Ø 3.0 mm X 20 cm
FT-204.48	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.91	Ø 2 mm X 15 cm
FT-206.92	Ø 2 mm X 20 cm
FT-206.93	Ø 2 mm X 25 cm
FT-206.94	Ø 2 mm X 30 cm
FT-206.95	Ø 2 mm X 35 cm
FT-206.96	Ø 3 mm X 15 cm
FT-206.97	Ø 3 mm X 20 cm
FT-206.98	Ø 3 mm X 25 cm
FT-206.99	Ø 3 mm X 30 cm
FT-206.100	Ø 3 mm X 35 cm
FT-206.101	Ø 4 mm X 15 cm
FT-206.102	Ø 4 mm X 20 cm
FT-206.103	Ø 4 mm X 25 cm
FT-206.104	Ø 4 mm X 30 cm
FT-206.105	Ø 4 mm X 35 cm
FT-206.106	Ø 5 mm X 15 cm
FT-206.107	Ø 5 mm X 20 cm
FT-206.108	Ø 5 mm X 25 cm
FT-206.109	Ø 5 mm X 30 cm
FT-206.110	Ø 5 mm X 35 cm
FT-206.111	Ø 6 mm X 15 cm
FT-206.112	Ø 6 mm X 20 cm
FT-206.113	Ø 6 mm X 25 cm
FT-206.114	Ø 6 mm X 30 cm
FT-206.115	Ø 6 mm X 35 cm



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula



Pyramid Holes Cannula

This commonly used design has three ports located 180 degrees from the thumb print. It is a slightly more aggressive cannula than the Standard during liposuction procedures.

60cc Tommeyhub

FT-201.41	Ø 2.5 mm X 15 cm
FT-201.42	Ø 2.5 mm X 20 cm
FT-201.43	Ø 3.0 mm X 15 cm
FT-201.44	Ø 3.0 mm X 20 cm
FT-201.45	Ø 3.0 mm X 25 cm
FT-201.46	Ø 4.0 mm X 25 cm
FT-201.47	Ø 4.0 mm X 30 cm
FT-201.48	Ø 4.0 mm X 35 cm
FT-201.49	Ø 4.5 mm X 30 cm
FT-201.50	Ø 4.5 mm X 35 cm

Treaded Hub

FT-202.41	Ø 3.0 mm X 15 cm
FT-202.42	Ø 3.0 mm X 25 cm
FT-202.43	Ø 4.0 mm X 15 cm
FT-202.44	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.41	Ø 1.5 mm X 12 cm
FT-203.42	Ø 2.0 mm X 15 cm
FT-203.43	Ø 2.5 mm X 15 cm
FT-203.44	Ø 2.5 mm X 20 cm
FT-203.45	Ø 3.0 mm X 15 cm
FT-203.46	Ø 3.0 mm X 20 cm
FT-203.47	Ø 3.0 mm X 25 cm

Luer lock

FT-204.51	Ø 1.5 mm X 12 cm
FT-204.52	Ø 2.0 mm X 15 cm
FT-204.53	Ø 2.5 mm X 15 cm
FT-204.54	Ø 2.5 mm X 20 cm
FT-204.55	Ø 3.0 mm X 15 cm
FT-204.56	Ø 3.0 mm X 20 cm
FT-204.57	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-205.51	Ø 1.5 mm X 12 cm
FT-205.52	Ø 2.0 mm X 15 cm
FT-205.53	Ø 2.0 mm X 20 cm
FT-205.54	Ø 2.5 mm X 15 cm
FT-205.55	Ø 2.5 mm X 20 cm
FT-205.56	Ø 3.0 mm X 15 cm
FT-205.57	Ø 3.0 mm X 20 cm
FT-205.58	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.121	Ø 2 mm X 15 cm
FT-206.122	Ø 2 mm X 20 cm
FT-206.123	Ø 2 mm X 25 cm
FT-206.124	Ø 2 mm X 30 cm
FT-206.125	Ø 2 mm X 35 cm
FT-206.126	Ø 3 mm X 15 cm
FT-206.127	Ø 3 mm X 20 cm
FT-206.128	Ø 3 mm X 25 cm
FT-206.129	Ø 3 mm X 30 cm
FT-206.130	Ø 3 mm X 35 cm
FT-206.131	Ø 4 mm X 15 cm
FT-206.132	Ø 4 mm X 20 cm
FT-206.133	Ø 4 mm X 25 cm
FT-206.134	Ø 4 mm X 30 cm
FT-206.135	Ø 4 mm X 35 cm
FT-206.136	Ø 5 mm X 15 cm
FT-206.137	Ø 5 mm X 20 cm
FT-206.138	Ø 5 mm X 25 cm
FT-206.139	Ø 5 mm X 30 cm
FT-206.140	Ø 5 mm X 35 cm
FT-206.141	Ø 6 mm X 15 cm
FT-206.142	Ø 6 mm X 20 cm
FT-206.143	Ø 6 mm X 25 cm
FT-206.144	Ø 6 mm X 30 cm
FT-206.145	Ø 6 mm X 35 cm



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula



Mercedes

This circumferential three port configuration cell friendly liposuction cannula is effective for multi-plane soft tissue harvesting.

60cc Tommeyhub

FT-201.51	Ø 2.5 mm X 15 cm
FT-201.52	Ø 2.5 mm X 20 cm
FT-201.53	Ø 3.0 mm X 15 cm
FT-201.54	Ø 3.0 mm X 20 cm
FT-201.55	Ø 3.0 mm X 25 cm
FT-201.56	Ø 4.0 mm X 25 cm
FT-201.57	Ø 4.0 mm X 30 cm
FT-201.58	Ø 4.0 mm X 35 cm
FT-201.59	Ø 4.5 mm X 30 cm
FT-201.60	Ø 4.5 mm X 35 cm

Luer lock (Stainless Steel)

FT-205.61	Ø 1.5 mm X 12 cm
FT-205.62	Ø 2.0 mm X 15 cm
FT-205.63	Ø 2.0 mm X 20 cm
FT-205.64	Ø 2.5 mm X 15 cm
FT-205.65	Ø 2.5 mm X 20 cm
FT-205.66	Ø 3.0 mm X 15 cm
FT-205.67	Ø 3.0 mm X 20 cm
FT-205.68	Ø 3.0 mm X 25 cm

Treaded Hub

FT-202.51	Ø 3.0 mm X 15 cm
FT-202.52	Ø 3.0 mm X 25 cm
FT-202.53	Ø 4.0 mm X 15 cm
FT-202.54	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.51	Ø 1.5 mm X 12 cm
FT-203.52	Ø 2.0 mm X 15 cm
FT-203.53	Ø 2.5 mm X 15 cm
FT-203.54	Ø 2.5 mm X 20 cm
FT-203.55	Ø 3.0 mm X 15 cm
FT-203.56	Ø 3.0 mm X 20 cm
FT-203.57	Ø 3.0 mm X 25 cm

Luer lock

FT-204.61	Ø 1.5 mm X 12 cm
FT-204.62	Ø 2.0 mm X 15 cm
FT-204.63	Ø 2.5 mm X 15 cm
FT-204.64	Ø 2.5 mm X 20 cm
FT-204.65	Ø 3.0 mm X 15 cm
FT-204.66	Ø 3.0 mm X 20 cm
FT-204.67	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.151	Ø 2 mm X 15 cm
FT-206.152	Ø 2 mm X 20 cm
FT-206.153	Ø 2 mm X 25 cm
FT-206.154	Ø 2 mm X 30 cm
FT-206.155	Ø 2 mm X 35 cm
FT-206.156	Ø 3 mm X 15 cm
FT-206.157	Ø 3 mm X 20 cm
FT-206.158	Ø 3 mm X 25 cm
FT-206.159	Ø 3 mm X 30 cm
FT-206.160	Ø 3 mm X 35 cm
FT-206.161	Ø 4 mm X 15 cm
FT-206.162	Ø 4 mm X 20 cm
FT-206.163	Ø 4 mm X 25 cm
FT-206.164	Ø 4 mm X 30 cm
FT-206.165	Ø 4 mm X 35 cm
FT-206.166	Ø 5 mm X 15 cm
FT-206.167	Ø 5 mm X 20 cm
FT-206.168	Ø 5 mm X 25 cm
FT-206.169	Ø 5 mm X 30 cm
FT-206.170	Ø 5 mm X 35 cm
FT-206.171	Ø 6 mm X 15 cm
FT-206.172	Ø 6 mm X 20 cm
FT-206.173	Ø 6 mm X 25 cm
FT-206.174	Ø 6 mm X 30 cm
FT-206.175	Ø 6 mm X 35 cm



3 Spiral Port Cannula (Carraway Style)

Spiral three port design on this Cell Friendly liposuction cannula is effective for multi-plane adipose tissue harvesting.

60ccTommeyhub

FT-201.61	Ø 2.5 mm X 15 cm
FT-201.62	Ø 2.5 mm X 20 cm
FT-201.63	Ø 3.0 mm X 15 cm
FT-201.64	Ø 3.0 mm X 20 cm
FT-201.65	Ø 3.0 mm X 25 cm
FT-201.66	Ø 4.0 mm X 25 cm
FT-201.67	Ø 4.0 mm X 30 cm
FT-201.68	Ø 4.0 mm X 35 cm
FT-201.69	Ø 4.5 mm X 30 cm
FT-201.70	Ø 4.5 mm X 35 cm

Treaded Hub

FT-202.61	Ø 3.0 mm X 15 cm
FT-202.62	Ø 3.0 mm X 25 cm
FT-202.63	Ø 4.0 mm X 15 cm
FT-202.64	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.61	Ø 1.5 mm X 12 cm
FT-203.62	Ø 2.0 mm X 15 cm
FT-203.63	Ø 2.5 mm X 15 cm
FT-203.64	Ø 2.5 mm X 20 cm
FT-203.65	Ø 3.0 mm X 15 cm
FT-203.66	Ø 3.0 mm X 20 cm
FT-203.67	Ø 3.0 mm X 25 cm

Luer lock


FT-204.71	Ø 1.5 mm X 12 cm
FT-204.72	Ø 2.0 mm X 15 cm
FT-204.73	Ø 2.5 mm X 15 cm
FT-204.74	Ø 2.5 mm X 20 cm
FT-204.75	Ø 3.0 mm X 15 cm
FT-204.76	Ø 3.0 mm X 20 cm
FT-204.77	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-205.71	Ø 1.5 mm X 12 cm
FT-205.72	Ø 2.0 mm X 15 cm
FT-205.73	Ø 2.0 mm X 20 cm
FT-205.74	Ø 2.5 mm X 15 cm
FT-205.75	Ø 2.5 mm X 20 cm
FT-205.76	Ø 3.0 mm X 15 cm
FT-205.77	Ø 3.0 mm X 20 cm
FT-205.78	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.181	Ø 2 mm X 15 cm
FT-206.182	Ø 2 mm X 20 cm
FT-206.183	Ø 2 mm X 25 cm
FT-206.184	Ø 2 mm X 30 cm
FT-206.185	Ø 2 mm X 35 cm
FT-206.186	Ø 3 mm X 15 cm
FT-206.187	Ø 3 mm X 20 cm
FT-206.188	Ø 3 mm X 25 cm
FT-206.189	Ø 3 mm X 30 cm
FT-206.190	Ø 3 mm X 35 cm
FT-206.191	Ø 4 mm X 15 cm
FT-206.192	Ø 4 mm X 20 cm
FT-206.193	Ø 4 mm X 25 cm
FT-206.194	Ø 4 mm X 30 cm
FT-206.195	Ø 4 mm X 35 cm
FT-206.196	Ø 5 mm X 15 cm
FT-206.197	Ø 5 mm X 20 cm
FT-206.198	Ø 5 mm X 25 cm
FT-206.199	Ø 5 mm X 30 cm
FT-206.200	Ø 5 mm X 35 cm
FT-206.201	Ø 6 mm X 15 cm
FT-206.202	Ø 6 mm X 20 cm
FT-206.203	Ø 6 mm X 25 cm
FT-206.204	Ø 6 mm X 30 cm
FT-206.205	Ø 6 mm X 35 cm



Six Spiral Port (Adjust below 3 Spiral Ports)

Six port spiral configuration on this Cell Friendly liposuction cannula harvests maximum tissue volume in minimal time.

60ccsmarthub

FT-201.91	Ø 2.5 mm X 15 cm
FT-201.92	Ø 2.5 mm X 20 cm
FT-201.93	Ø 3.0 mm X 15 cm
FT-201.94	Ø 3.0 mm X 25 cm
FT-201.95	Ø 4.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-205.231	Ø 2.5 mm X 15 cm
FT-205.232	Ø 3.0 mm X 15 cm
FT-205.233	Ø 3.0 mm X 25 cm
FT-205.234	Ø 4.0 mm X 15 cm
FT-205.235	Ø 4.0 mm X 25 cm

Treaded Hub

FT-202.81	Ø 3.0 mm X 15 cm
FT-202.82	Ø 3.0 mm X 25 cm
FT-202.83	Ø 4.0 mm X 15 cm
FT-202.84	Ø 4.0 mm X 25 cm

Power Cannula

FT-206.901	2mm OD X 15 cm
FT-206.902	2mm OD X 20 cm
FT-206.903	2mm OD X 25 cm
FT-206.904	2mm OD X 30 cm
FT-206.905	2mm OD X 35 cm
FT-206.906	3mm OD X 15 cm
FT-206.907	3mm OD X 20 cm
FT-206.908	3mm OD X 25 cm
FT-206.909	3mm OD X 30 cm
FT-206.910	3mm OD X 35 cm
FT-206.911	4mm OD X 15 cm
FT-206.912	4mm OD X 20 cm
FT-206.913	4mm OD X 25 cm
FT-206.914	4mm OD X 30 cm
FT-206.915	4mm OD X 35 cm
FT-206.916	5mm OD X 15 cm
FT-206.917	5mm OD X 20 cm
FT-206.918	5mm OD X 25 cm
FT-206.919	5mm OD X 30 cm
FT-206.920	5mm OD X 35 cm
FT-206.921	6mm OD X 15 cm
FT-206.922	6mm OD X 20 cm
FT-206.923	6mm OD X 25 cm
FT-206.924	6mm OD X 30 cm
FT-206.925	6mm OD X 35 cm

Ultra Luer lock

FT-203.81	Ø 2.5 mm X 15 cm
FT-203.82	Ø 3.0 mm X 15 cm
FT-203.83	Ø 3.0 mm X 25 cm
FT-203.84	Ø 4.0 mm X 15 cm
FT-203.85	Ø 4.0 mm X 25 cm

Luer lock

FT-204.91	Ø 2.5 mm X 15 cm
FT-204.92	Ø 3.0 mm X 15 cm
FT-204.93	Ø 3.0 mm X 25 cm
FT-204.94	Ø 4.0 mm X 15 cm
FT-204.95	Ø 4.0 mm X 25 cm



Luer lock

Luer lock
(Stainless Steel)

Ultra Luer lock

Treaded

60cc Tommey hub

Power Cannula



Three Port Cannula

This body sculpting instrument is excellent for harvesting small aliquots of tissue during liposuction procedures.

60cc Tommeyhub

FT-201.71	Ø 2.5 mm X 15 cm
FT-201.72	Ø 2.5 mm X 20 cm
FT-201.73	Ø 3.0 mm X 15 cm
FT-201.74	Ø 3.0 mm X 20 cm
FT-201.75	Ø 3.0 mm X 25 cm
FT-201.76	Ø 4.0 mm X 25 cm
FT-201.77	Ø 4.0 mm X 30 cm
FT-201.78	Ø 4.0 mm X 35 cm
FT-201.79	Ø 4.5 mm X 30 cm
FT-201.80	Ø 4.5 mm X 35 cm

Treaded Hub

FT-202.71	Ø 3.0 mm X 15 cm
FT-202.72	Ø 3.0 mm X 25 cm
FT-202.73	Ø 4.0 mm X 15 cm
FT-202.74	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-203.71	Ø 1.5 mm X 12 cm
FT-203.72	Ø 2.0 mm X 15 cm
FT-203.73	Ø 2.5 mm X 15 cm
FT-203.74	Ø 2.5 mm X 20 cm
FT-203.75	Ø 3.0 mm X 15 cm
FT-203.76	Ø 3.0 mm X 20 cm
FT-203.77	Ø 3.0 mm X 25 cm

Luer lock


FT-204.81	Ø 1.5 mm X 12 cm
FT-204.82	Ø 2.0 mm X 15 cm
FT-204.83	Ø 2.5 mm X 15 cm
FT-204.84	Ø 2.5 mm X 20 cm
FT-204.85	Ø 3.0 mm X 15 cm
FT-204.86	Ø 3.0 mm X 20 cm
FT-204.87	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-205.81	Ø 1.5 mm X 12 cm
FT-205.82	Ø 2.0 mm X 15 cm
FT-205.83	Ø 2.0 mm X 20 cm
FT-205.84	Ø 2.5 mm X 15 cm
FT-205.85	Ø 2.5 mm X 20 cm
FT-205.86	Ø 3.0 mm X 15 cm
FT-205.87	Ø 3.0 mm X 20 cm
FT-205.88	Ø 3.0 mm X 25 cm

Power Cannula

FT-206.211	Ø 2 mm X 15 cm
FT-206.212	Ø 2 mm X 20 cm
FT-206.213	Ø 2 mm X 25 cm
FT-206.214	Ø 2 mm X 30 cm
FT-206.215	Ø 2 mm X 35 cm
FT-206.216	Ø 3 mm X 15 cm
FT-206.217	Ø 3 mm X 20 cm
FT-206.218	Ø 3 mm X 25 cm
FT-206.219	Ø 3 mm X 30 cm
FT-206.220	Ø 3 mm X 35 cm
FT-206.221	Ø 4 mm X 15 cm
FT-206.222	Ø 4 mm X 20 cm
FT-206.223	Ø 4 mm X 25 cm
FT-206.224	Ø 4 mm X 30 cm
FT-206.225	Ø 4 mm X 35 cm
FT-206.226	Ø 5 mm X 15 cm
FT-206.227	Ø 5 mm X 20 cm
FT-206.228	Ø 5 mm X 25 cm
FT-206.229	Ø 5 mm X 30 cm
FT-206.230	Ø 5 mm X 35 cm
FT-206.231	Ø 6 mm X 15 cm
FT-206.232	Ø 6 mm X 20 cm
FT-206.233	Ø 6 mm X 25 cm
FT-206.234	Ø 6 mm X 30 cm
FT-206.235	Ø 6 mm X 35 cm



Basket Shaped Cannula

The Cell Friendly Basket Cannula is a versatile liposculpting tool. It can operate as an equalizer without vacuum to treat adipose tissue and skin irregularities in revision liposuction and fat transfer procedures. The Basket Cannula is also adept under vacuum as a fat harvesting cannula. For finishing, this cannula can be used as an injector for equalizing and adding adipose tissue to sculpt smoother contours.

60cc Tommeyhub

FT-201.81	Ø 3.0 mm X 15 cm
FT-201.82	Ø 3.0 mm X 20 cm
FT-201.83	Ø 3.0 mm X 25 cm
FT-201.84	Ø 4.0 mm X 25 cm
FT-201.85	Ø 4.0 mm X 30 cm
FT-201.86	Ø 4.0 mm X 35 cm



Hole on Tip-Style I

Power Cannula

FT-206.251	Ø 2 mm X 15 cm	FT-206.260	Ø 3 mm X 35 cm	FT-206.268	Ø 5 mm X 25 cm
FT-206.252	Ø 2 mm X 20 cm	FT-206.261	Ø 4 mm X 15 cm	FT-206.269	Ø 5 mm X 30 cm
FT-206.253	Ø 2 mm X 25 cm	FT-206.262	Ø 4 mm X 20 cm	FT-206.270	Ø 5 mm X 35 cm
FT-206.254	Ø 2 mm X 30 cm	FT-206.263	Ø 4 mm X 25 cm	FT-206.271	Ø 6 mm X 15 cm
FT-206.255	Ø 2 mm X 35 cm	FT-206.264	Ø 4 mm X 30 cm	FT-206.272	Ø 6 mm X 20 cm
FT-206.256	Ø 3 mm X 15 cm	FT-206.265	Ø 4 mm X 35 cm	FT-206.273	Ø 6 mm X 25 cm
FT-206.257	Ø 3 mm X 20 cm	FT-206.266	Ø 5 mm X 15 cm	FT-206.274	Ø 6 mm X 30 cm
FT-206.258	Ø 3 mm X 25 cm	FT-206.267	Ø 5 mm X 20 cm	FT-206.275	Ø 6 mm X 35 cm
FT-206.259	Ø 3 mm X 30 cm				



Holeontip-Style II

Power Cannula

FT-206.281	Ø 2 mm X 15 cm	FT-206.290	Ø 3 mm X 35 cm	FT-206.298	Ø 5 mm X 25 cm
FT-206.282	Ø 2 mm X 20 cm	FT-206.291	Ø 4 mm X 15 cm	FT-206.299	Ø 5 mm X 30 cm
FT-206.283	Ø 2 mm X 25 cm	FT-206.292	Ø 4 mm X 20 cm	FT-206.300	Ø 5 mm X 35 cm
FT-206.284	Ø 2 mm X 30 cm	FT-206.293	Ø 4 mm X 25 cm	FT-206.301	Ø 6 mm X 15 cm
FT-206.285	Ø 2 mm X 35 cm	FT-206.294	Ø 4 mm X 30 cm	FT-206.302	Ø 6 mm X 20 cm
FT-206.286	Ø 3 mm X 15 cm	FT-206.295	Ø 4 mm X 35 cm	FT-206.303	Ø 6 mm X 25 cm
FT-206.287	Ø 3 mm X 20 cm	FT-206.296	Ø 5 mm X 15 cm	FT-206.304	Ø 6 mm X 30 cm
FT-206.288	Ø 3 mm X 25 cm	FT-206.297	Ø 5 mm X 20 cm	FT-206.305	Ø 6 mm X 35 cm
FT-206.289	Ø 3 mm X 30 cm				



Deep Sharp Cut Dissector

Power Cannula

FT-206.311	Ø 2 mm X 15 cm	FT-206.320	Ø 3 mm X 35 cm	FT-206.328	Ø 5 mm X 25 cm
FT-206.312	Ø 2 mm X 20 cm	FT-206.321	Ø 4 mm X 15 cm	FT-206.329	Ø 5 mm X 30 cm
FT-206.313	Ø 2 mm X 25 cm	FT-206.322	Ø 4 mm X 20 cm	FT-206.330	Ø 5 mm X 35 cm
FT-206.314	Ø 2 mm X 30 cm	FT-206.323	Ø 4 mm X 25 cm	FT-206.331	Ø 6 mm X 15 cm
FT-206.315	Ø 2 mm X 35 cm	FT-206.324	Ø 4 mm X 30 cm	FT-206.332	Ø 6 mm X 20 cm
FT-206.316	Ø 3 mm X 15 cm	FT-206.325	Ø 4 mm X 35 cm	FT-206.333	Ø 6 mm X 25 cm
FT-206.317	Ø 3 mm X 20 cm	FT-206.326	Ø 5 mm X 15 cm	FT-206.334	Ø 6 mm X 30 cm
FT-206.318	Ø 3 mm X 25 cm	FT-206.327	Ø 5 mm X 20 cm	FT-206.335	Ø 6 mm X 35 cm
FT-206.319	Ø 3 mm X 30 cm				



Bi-Port Cannula

Luer lock (Stainless Steel)

	10cm	15cm	20cm	25cm
2.0 Ø	FT-205.91	FT-205.101	FT-205.111	FT-205.121
2.5 Ø	FT-205.92	FT-205.102	FT-205.112	FT-205.122
3.0 Ø	FT-205.93	FT-205.103	FT-205.113	FT-205.123
3.5 Ø	FT-205.94	FT-205.104	FT-205.114	FT-205.124
4.0 Ø	FT-205.95	FT-205.105	FT-205.115	FT-205.125

Power Cannula

FT-206.351	Ø 2 mm X 15 cm	FT-206.360	Ø 3 mm X 35 cm	FT-206.368	Ø 5 mm X 25 cm
FT-206.352	Ø 2 mm X 20 cm	FT-206.361	Ø 4 mm X 15 cm	FT-206.369	Ø 5 mm X 30 cm
FT-206.353	Ø 2 mm X 25 cm	FT-206.362	Ø 4 mm X 20 cm	FT-206.370	Ø 5 mm X 35 cm
FT-206.354	Ø 2 mm X 30 cm	FT-206.363	Ø 4 mm X 25 cm	FT-206.371	Ø 6 mm X 15 cm
FT-206.355	Ø 2 mm X 35 cm	FT-206.364	Ø 4 mm X 30 cm	FT-206.372	Ø 6 mm X 20 cm
FT-206.356	Ø 3 mm X 15 cm	FT-206.365	Ø 4 mm X 35 cm	FT-206.373	Ø 6 mm X 25 cm
FT-206.357	Ø 3 mm X 20 cm	FT-206.366	Ø 5 mm X 15 cm	FT-206.374	Ø 6 mm X 30 cm
FT-206.358	Ø 3 mm X 25 cm	FT-206.367	Ø 5 mm X 20 cm	FT-206.375	Ø 6 mm X 35 cm
FT-206.359	Ø 3 mm X 30 cm				



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula



Tapered Tip Single Central Hole

Power Cannula

FT-206.391	Ø 2 mm X 15 cm	FT-206.400	Ø 3 mm X 35 cm	FT-206.408	Ø 5 mm X 25 cm
FT-206.392	Ø 2 mm X 20 cm	FT-206.401	Ø 4 mm X 15 cm	FT-206.409	Ø 5 mm X 30 cm
FT-206.393	Ø 2 mm X 25 cm	FT-206.402	Ø 4 mm X 20 cm	FT-206.410	Ø 5 mm X 35 cm
FT-206.394	Ø 2 mm X 30 cm	FT-206.403	Ø 4 mm X 25 cm	FT-206.411	Ø 6 mm X 15 cm
FT-206.395	Ø 2 mm X 35 cm	FT-206.404	Ø 4 mm X 30 cm	FT-206.412	Ø 6 mm X 20 cm
FT-206.396	Ø 3 mm X 15 cm	FT-206.405	Ø 4 mm X 35 cm	FT-206.413	Ø 6 mm X 25 cm
FT-206.397	Ø 3 mm X 20 cm	FT-206.406	Ø 5 mm X 15 cm	FT-206.414	Ø 6 mm X 30 cm
FT-206.398	Ø 3 mm X 25 cm	FT-206.407	Ø 5 mm X 20 cm	FT-206.415	Ø 6 mm X 35 cm
FT-206.399	Ø 3 mm X 30 cm				



Taperd Tip Two Spiral Holes

Power Cannula

FT-206.431	Ø 2 mm X 15 cm	FT-206.440	Ø 3 mm X 35 cm	FT-206.448	Ø 5 mm X 25 cm
FT-206.432	Ø 2 mm X 20 cm	FT-206.441	Ø 4 mm X 15 cm	FT-206.449	Ø 5 mm X 30 cm
FT-206.433	Ø 2 mm X 25 cm	FT-206.442	Ø 4 mm X 20 cm	FT-206.450	Ø 5 mm X 35 cm
FT-206.434	Ø 2 mm X 30 cm	FT-206.443	Ø 4 mm X 25 cm	FT-206.451	Ø 6 mm X 15 cm
FT-206.435	Ø 2 mm X 35 cm	FT-206.444	Ø 4 mm X 30 cm	FT-206.452	Ø 6 mm X 20 cm
FT-206.436	Ø 3 mm X 15 cm	FT-206.445	Ø 4 mm X 35 cm	FT-206.453	Ø 6 mm X 25 cm
FT-206.437	Ø 3 mm X 20 cm	FT-206.446	Ø 5 mm X 15 cm	FT-206.454	Ø 6 mm X 30 cm
FT-206.438	Ø 3 mm X 25 cm	FT-206.447	Ø 5 mm X 20 cm	FT-206.455	Ø 6 mm X 35 cm
FT-206.439	Ø 3 mm X 30 cm				



Fig. 3F

Power Cannula

FT-206.471	Ø 2 mm X 15 cm	FT-206.480	Ø 3 mm X 35 cm	FT-206.488	Ø 5 mm X 25 cm
FT-206.472	Ø 2 mm X 20 cm	FT-206.481	Ø 4 mm X 15 cm	FT-206.489	Ø 5 mm X 30 cm
FT-206.473	Ø 2 mm X 25 cm	FT-206.482	Ø 4 mm X 20 cm	FT-206.490	Ø 5 mm X 35 cm
FT-206.474	Ø 2 mm X 30 cm	FT-206.483	Ø 4 mm X 25 cm	FT-206.491	Ø 6 mm X 15 cm
FT-206.475	Ø 2 mm X 35 cm	FT-206.484	Ø 4 mm X 30 cm	FT-206.492	Ø 6 mm X 20 cm
FT-206.476	Ø 3 mm X 15 cm	FT-206.485	Ø 4 mm X 35 cm	FT-206.493	Ø 6 mm X 25 cm
FT-206.477	Ø 3 mm X 20 cm	FT-206.486	Ø 5 mm X 15 cm	FT-206.494	Ø 6 mm X 30 cm
FT-206.478	Ø 3 mm X 25 cm	FT-206.487	Ø 5 mm X 20 cm	FT-206.495	Ø 6 mm X 35 cm




Fig.2Ech

Power Cannula

FT-206.501	Ø 2 mm X 15 cm	FT-206.510	Ø 3 mm X 35 cm	FT-206.518	Ø 5 mm X 25 cm
FT-206.502	Ø 2 mm X 20 cm	FT-206.511	Ø 4 mm X 15 cm	FT-206.519	Ø 5 mm X 30 cm
FT-206.503	Ø 2 mm X 25 cm	FT-206.512	Ø 4 mm X 20 cm	FT-206.520	Ø 5 mm X 35 cm
FT-206.504	Ø 2 mm X 30 cm	FT-206.513	Ø 4 mm X 25 cm	FT-206.521	Ø 6 mm X 15 cm
FT-206.505	Ø 2 mm X 35 cm	FT-206.514	Ø 4 mm X 30 cm	FT-206.522	Ø 6 mm X 20 cm
FT-206.506	Ø 3 mm X 15 cm	FT-206.515	Ø 4 mm X 35 cm	FT-206.523	Ø 6 mm X 25 cm
FT-206.507	Ø 3 mm X 20 cm	FT-206.516	Ø 5 mm X 15 cm	FT-206.524	Ø 6 mm X 30 cm
FT-206.508	Ø 3 mm X 25 cm	FT-206.517	Ø 5 mm X 20 cm	FT-206.525	Ø 6 mm X 35 cm
FT-206.509	Ø 3 mm X 30 cm				



Two Square-In Line

Luer lock (Stainless Steel)

	7cm	10cm	15cm	20cm	25cm
2.0 Ø	FT-205.131	FT-205.141	FT-205.151	FT-205.161	FT-205.171
2.5 Ø	FT-205.132	FT-205.142	FT-205.152	FT-205.162	FT-205.172
3.0 Ø	FT-205.133	FT-205.143	FT-205.153	FT-205.163	FT-205.173
3.5 Ø	FT-205.134	FT-205.144	FT-205.154	FT-205.164	FT-205.174
4.0 Ø	FT-205.135	FT-205.145	FT-205.155	FT-205.165	FT-205.175

Power Cannula

FT-206.551	Ø 2 mm X 15 cm	FT-206.560	Ø 3 mm X 35 cm	FT-206.568	Ø 5 mm X 25 cm
FT-206.552	Ø 2 mm X 20 cm	FT-206.561	Ø 4 mm X 15 cm	FT-206.569	Ø 5 mm X 30 cm
FT-206.553	Ø 2 mm X 25 cm	FT-206.562	Ø 4 mm X 20 cm	FT-206.570	Ø 5 mm X 35 cm
FT-206.554	Ø 2 mm X 30 cm	FT-206.563	Ø 4 mm X 25 cm	FT-206.571	Ø 6 mm X 15 cm
FT-206.555	Ø 2 mm X 35 cm	FT-206.564	Ø 4 mm X 30 cm	FT-206.572	Ø 6 mm X 20 cm
FT-206.556	Ø 3 mm X 15 cm	FT-206.565	Ø 4 mm X 35 cm	FT-206.573	Ø 6 mm X 25 cm
FT-206.557	Ø 3 mm X 20 cm	FT-206.566	Ø 5 mm X 15 cm	FT-206.574	Ø 6 mm X 30 cm
FT-206.558	Ø 3 mm X 25 cm	FT-206.567	Ø 5 mm X 20 cm	FT-206.575	Ø 6 mm X 35 cm
FT-206.559	Ø 3 mm X 30 cm				



One Central Port-Two Parallel Ports

Power Cannula

FT-206.591	Ø 2 mm X 15 cm	FT-206.600	Ø 3 mm X 35 cm	FT-206.608	Ø 5 mm X 25 cm
FT-206.592	Ø 2 mm X 20 cm	FT-206.601	Ø 4 mm X 15 cm	FT-206.609	Ø 5 mm X 30 cm
FT-206.593	Ø 2 mm X 25 cm	FT-206.602	Ø 4 mm X 20 cm	FT-206.610	Ø 5 mm X 35 cm
FT-206.594	Ø 2 mm X 30 cm	FT-206.603	Ø 4 mm X 25 cm	FT-206.611	Ø 6 mm X 15 cm
FT-206.595	Ø 2 mm X 35 cm	FT-206.604	Ø 4 mm X 30 cm	FT-206.612	Ø 6 mm X 20 cm
FT-206.596	Ø 3 mm X 15 cm	FT-206.605	Ø 4 mm X 35 cm	FT-206.613	Ø 6 mm X 25 cm
FT-206.597	Ø 3 mm X 20 cm	FT-206.606	Ø 5 mm X 15 cm	FT-206.614	Ø 6 mm X 30 cm
FT-206.598	Ø 3 mm X 25 cm	FT-206.607	Ø 5 mm X 20 cm	FT-206.615	Ø 6 mm X 35 cm
FT-206.599	Ø 3 mm X 30 cm				



Two Parallel Port- One Central Port

Power Cannula

FT-206.651	Ø 2 mm X 15 cm	FT-206.660	Ø 3 mm X 35 cm	FT-206.668	Ø 5 mm X 25 cm
FT-206.652	Ø 2 mm X 20 cm	FT-206.661	Ø 4 mm X 15 cm	FT-206.669	Ø 5 mm X 30 cm
FT-206.653	Ø 2 mm X 25 cm	FT-206.662	Ø 4 mm X 20 cm	FT-206.670	Ø 5 mm X 35 cm
FT-206.654	Ø 2 mm X 30 cm	FT-206.663	Ø 4 mm X 25 cm	FT-206.671	Ø 6 mm X 15 cm
FT-206.655	Ø 2 mm X 35 cm	FT-206.664	Ø 4 mm X 30 cm	FT-206.672	Ø 6 mm X 20 cm
FT-206.656	Ø 3 mm X 15 cm	FT-206.665	Ø 4 mm X 35 cm	FT-206.673	Ø 6 mm X 25 cm
FT-206.657	Ø 3 mm X 20 cm	FT-206.666	Ø 5 mm X 15 cm	FT-206.674	Ø 6 mm X 30 cm
FT-206.658	Ø 3 mm X 25 cm	FT-206.667	Ø 5 mm X 20 cm	FT-206.675	Ø 6 mm X 35 cm
FT-206.659	Ø 3 mm X 30 cm				



Fig. 3Me

Luer lock (Stainless Steel)

	7cm	10cm	15cm	20cm	25cm
2.0 Ø	FT-205.181	FT-205.191	FT-205.201	FT-205.211	FT-205.221
2.5 Ø	FT-205.182	FT-205.192	FT-205.202	FT-205.212	FT-205.222
3.0 Ø	FT-205.183	FT-205.193	FT-205.203	FT-205.213	FT-205.223
3.5 Ø	FT-205.184	FT-205.194	FT-205.204	FT-205.214	FT-205.224
4.0 Ø	FT-205.185	FT-205.195	FT-205.205	FT-205.215	FT-205.225

Power Cannula

FT-206.681	Ø 2 mm X 15 cm	FT-206.690	Ø 3 mm X 35 cm	FT-206.698	Ø 5 mm X 25 cm
FT-206.682	Ø 2 mm X 20 cm	FT-206.691	Ø 4 mm X 15 cm	FT-206.699	Ø 5 mm X 30 cm
FT-206.683	Ø 2 mm X 25 cm	FT-206.692	Ø 4 mm X 20 cm	FT-206.700	Ø 5 mm X 35 cm
FT-206.684	Ø 2 mm X 30 cm	FT-206.693	Ø 4 mm X 25 cm	FT-206.701	Ø 6 mm X 15 cm
FT-206.685	Ø 2 mm X 35 cm	FT-206.694	Ø 4 mm X 30 cm	FT-206.702	Ø 6 mm X 20 cm
FT-206.686	Ø 3 mm X 15 cm	FT-206.695	Ø 4 mm X 35 cm	FT-206.703	Ø 6 mm X 25 cm
FT-206.687	Ø 3 mm X 20 cm	FT-206.696	Ø 5 mm X 15 cm	FT-206.704	Ø 6 mm X 30 cm
FT-206.688	Ø 3 mm X 25 cm	FT-206.697	Ø 5 mm X 20 cm	FT-206.705	Ø 6 mm X 35 cm
FT-206.689	Ø 3 mm X 30 cm				



Fig. 3Hr

Power Cannula

FT-206.721	Ø 2 mm X 15 cm	FT-206.730	Ø 3 mm X 35 cm	FT-206.738	Ø 5 mm X 25 cm
FT-206.722	Ø 2 mm X 20 cm	FT-206.731	Ø 4 mm X 15 cm	FT-206.739	Ø 5 mm X 30 cm
FT-206.723	Ø 2 mm X 25 cm	FT-206.732	Ø 4 mm X 20 cm	FT-206.740	Ø 5 mm X 35 cm
FT-206.724	Ø 2 mm X 30 cm	FT-206.733	Ø 4 mm X 25 cm	FT-206.741	Ø 6 mm X 15 cm
FT-206.725	Ø 2 mm X 35 cm	FT-206.734	Ø 4 mm X 30 cm	FT-206.742	Ø 6 mm X 20 cm
FT-206.726	Ø 3 mm X 15 cm	FT-206.735	Ø 4 mm X 35 cm	FT-206.743	Ø 6 mm X 25 cm
FT-206.727	Ø 3 mm X 20 cm	FT-206.736	Ø 5 mm X 15 cm	FT-206.744	Ø 6 mm X 30 cm
FT-206.728	Ø 3 mm X 25 cm	FT-206.737	Ø 5 mm X 20 cm	FT-206.745	Ø 6 mm X 35 cm
FT-206.729	Ø 3 mm X 30 cm				



Fig.Riz

Power Cannula

FT-206.751	Ø 2 mm X 15 cm	FT-206.760	Ø 3 mm X 35 cm	FT-206.768	Ø 5 mm X 25 cm
FT-206.752	Ø 2 mm X 20 cm	FT-206.761	Ø 4 mm X 15 cm	FT-206.769	Ø 5 mm X 30 cm
FT-206.753	Ø 2 mm X 25 cm	FT-206.762	Ø 4 mm X 20 cm	FT-206.770	Ø 5 mm X 35 cm
FT-206.754	Ø 2 mm X 30 cm	FT-206.763	Ø 4 mm X 25 cm	FT-206.771	Ø 6 mm X 15 cm
FT-206.755	Ø 2 mm X 35 cm	FT-206.764	Ø 4 mm X 30 cm	FT-206.772	Ø 6 mm X 20 cm
FT-206.756	Ø 3 mm X 15 cm	FT-206.765	Ø 4 mm X 35 cm	FT-206.773	Ø 6 mm X 25 cm
FT-206.757	Ø 3 mm X 20 cm	FT-206.766	Ø 5 mm X 15 cm	FT-206.774	Ø 6 mm X 30 cm
FT-206.758	Ø 3 mm X 25 cm	FT-206.767	Ø 5 mm X 20 cm	FT-206.775	Ø 6 mm X 35 cm
FT-206.759	Ø 3 mm X 30 cm				




Fig. 5Ai

Power Cannula

FT-206.721	Ø 2 mm X 15 cm	FT-206.730	Ø 3 mm X 35 cm	FT-206.738	Ø 5 mm X 25 cm
FT-206.722	Ø 2 mm X 20 cm	FT-206.731	Ø 4 mm X 15 cm	FT-206.739	Ø 5 mm X 30 cm
FT-206.723	Ø 2 mm X 25 cm	FT-206.732	Ø 4 mm X 20 cm	FT-206.740	Ø 5 mm X 35 cm
FT-206.724	Ø 2 mm X 30 cm	FT-206.733	Ø 4 mm X 25 cm	FT-206.741	Ø 6 mm X 15 cm
FT-206.725	Ø 2 mm X 35 cm	FT-206.734	Ø 4 mm X 30 cm	FT-206.742	Ø 6 mm X 20 cm
FT-206.726	Ø 3 mm X 15 cm	FT-206.735	Ø 4 mm X 35 cm	FT-206.743	Ø 6 mm X 25 cm
FT-206.727	Ø 3 mm X 20 cm	FT-206.736	Ø 5 mm X 15 cm	FT-206.744	Ø 6 mm X 30 cm
FT-206.728	Ø 3 mm X 25 cm	FT-206.737	Ø 5 mm X 20 cm	FT-206.745	Ø 6 mm X 35 cm
FT-206.729	Ø 3 mm X 30 cm				



Fig. Mr

Power Cannula

FT-206.751	Ø 2 mm X 15 cm	FT-206.760	Ø 3 mm X 35 cm	FT-206.768	Ø 5 mm X 25 cm
FT-206.752	Ø 2 mm X 20 cm	FT-206.761	Ø 4 mm X 15 cm	FT-206.769	Ø 5 mm X 30 cm
FT-206.753	Ø 2 mm X 25 cm	FT-206.762	Ø 4 mm X 20 cm	FT-206.770	Ø 5 mm X 35 cm
FT-206.754	Ø 2 mm X 30 cm	FT-206.763	Ø 4 mm X 25 cm	FT-206.771	Ø 6 mm X 15 cm
FT-206.755	Ø 2 mm X 35 cm	FT-206.764	Ø 4 mm X 30 cm	FT-206.772	Ø 6 mm X 20 cm
FT-206.756	Ø 3 mm X 15 cm	FT-206.765	Ø 4 mm X 35 cm	FT-206.773	Ø 6 mm X 25 cm
FT-206.757	Ø 3 mm X 20 cm	FT-206.766	Ø 5 mm X 15 cm	FT-206.774	Ø 6 mm X 30 cm
FT-206.758	Ø 3 mm X 25 cm	FT-206.767	Ø 5 mm X 20 cm	FT-206.775	Ø 6 mm X 35 cm
FT-206.759	Ø 3 mm X 30 cm				



Fig. P3G

Power Cannula

FT-206.791	Ø 2 mm X 15 cm	FT-206.800	Ø 3 mm X 35 cm	FT-206.808	Ø 5 mm X 25 cm
FT-206.792	Ø 2 mm X 20 cm	FT-206.801	Ø 4 mm X 15 cm	FT-206.809	Ø 5 mm X 30 cm
FT-206.793	Ø 2 mm X 25 cm	FT-206.802	Ø 4 mm X 20 cm	FT-206.810	Ø 5 mm X 35 cm
FT-206.794	Ø 2 mm X 30 cm	FT-206.803	Ø 4 mm X 25 cm	FT-206.811	Ø 6 mm X 15 cm
FT-206.795	Ø 2 mm X 35 cm	FT-206.804	Ø 4 mm X 30 cm	FT-206.812	Ø 6 mm X 20 cm
FT-206.796	Ø 3 mm X 15 cm	FT-206.805	Ø 4 mm X 35 cm	FT-206.813	Ø 6 mm X 25 cm
FT-206.797	Ø 3 mm X 20 cm	FT-206.806	Ø 5 mm X 15 cm	FT-206.814	Ø 6 mm X 30 cm
FT-206.798	Ø 3 mm X 25 cm	FT-206.807	Ø 5 mm X 20 cm	FT-206.815	Ø 6 mm X 35 cm
FT-206.799	Ø 3 mm X 30 cm				




Fig.Py

Power Cannula					
FT-206.831	Ø 2 mm X 15 cm	FT-206.840	Ø 3 mm X 35 cm	FT-206.848	Ø 5 mm X 25 cm
FT-206.832	Ø 2 mm X 20 cm	FT-206.841	Ø 4 mm X 15 cm	FT-206.849	Ø 5 mm X 30 cm
FT-206.833	Ø 2 mm X 25 cm	FT-206.842	Ø 4 mm X 20 cm	FT-206.850	Ø 5 mm X 35 cm
FT-206.834	Ø 2 mm X 30 cm	FT-206.843	Ø 4 mm X 25 cm	FT-206.851	Ø 6 mm X 15 cm
FT-206.835	Ø 2 mm X 35 cm	FT-206.844	Ø 4 mm X 30 cm	FT-206.852	Ø 6 mm X 20 cm
FT-206.836	Ø 3 mm X 15 cm	FT-206.845	Ø 4 mm X 35 cm	FT-206.853	Ø 6 mm X 25 cm
FT-206.837	Ø 3 mm X 20 cm	FT-206.846	Ø 5 mm X 15 cm	FT-206.854	Ø 6 mm X 30 cm
FT-206.838	Ø 3 mm X 25 cm	FT-206.847	Ø 5 mm X 20 cm	FT-206.855	Ø 6 mm X 35 cm
FT-206.839	Ø 3 mm X 30 cm				



Fig.Le

Power Cannula					
FT-206.861	Ø 2 mm X 15 cm	FT-206.870	Ø 3 mm X 35 cm	FT-206.878	Ø 5 mm X 25 cm
FT-206.862	Ø 2 mm X 20 cm	FT-206.871	Ø 4 mm X 15 cm	FT-206.879	Ø 5 mm X 30 cm
FT-206.863	Ø 2 mm X 25 cm	FT-206.872	Ø 4 mm X 20 cm	FT-206.880	Ø 5 mm X 35 cm
FT-206.864	Ø 2 mm X 30 cm	FT-206.873	Ø 4 mm X 25 cm	FT-206.881	Ø 6 mm X 15 cm
FT-206.865	Ø 2 mm X 35 cm	FT-206.874	Ø 4 mm X 30 cm	FT-206.882	Ø 6 mm X 20 cm
FT-206.866	Ø 3 mm X 15 cm	FT-206.875	Ø 4 mm X 35 cm	FT-206.883	Ø 6 mm X 25 cm
FT-206.867	Ø 3 mm X 20 cm	FT-206.876	Ø 5 mm X 15 cm	FT-206.884	Ø 6 mm X 30 cm
FT-206.868	Ø 3 mm X 25 cm	FT-206.877	Ø 5 mm X 20 cm	FT-206.885	Ø 6 mm X 35 cm
FT-206.869	Ø 3 mm X 30 cm				



Harvester Cannulas

Fat Harvesting for Re-grafting is a delicate difficult and time consuming process. But with HELIX Surgical's Specifically designed Harvesting Cannulas will help you to achieve this goal in super efficient way.

Micro fat Harvesting - No problem with HELIX Cannulas.



Multiport Speed Harvester- Both End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

60ccTommeyhub

FT-301.01	Ø 2.5 mm X 15 cm
FT-301.02	Ø 2.5 mm X 20 cm
FT-301.03	Ø 3.0 mm X 15 cm
FT-301.04	Ø 3.0 mm X 20 cm
FT-301.05	Ø 3.0 mm X 25 cm
FT-301.06	Ø 4.0 mm X 25 cm
FT-301.07	Ø 4.0 mm X 30 cm
FT-301.08	Ø 4.0 mm X 35 cm
FT-301.09	Ø 4.5 mm X 30 cm
FT-301.10	Ø 4.5 mm X 35 cm

Treaded Hub

FT-302.01	Ø 3.0 mm X 15 cm
FT-302.02	Ø 3.0 mm X 25 cm
FT-302.03	Ø 4.0 mm X 15 cm
FT-302.04	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-303.01	Ø 2.0 mm X 15 cm
FT-303.02	Ø 2.5 mm X 15 cm
FT-303.03	Ø 2.5 mm X 20 cm
FT-303.04	Ø 3.0 mm X 15 cm
FT-303.05	Ø 3.0 mm X 20 cm
FT-303.06	Ø 3.0 mm X 25 cm

Luer lock

FT-304.01	Ø 2.0 mm X 15 cm
FT-304.02	Ø 2.5 mm X 15 cm
FT-304.03	Ø 2.5 mm X 20 cm
FT-304.04	Ø 3.0 mm X 15 cm
FT-304.05	Ø 3.0 mm X 20 cm
FT-304.06	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-305.01	Ø 2.0 mm X 15 cm
FT-305.02	Ø 2.5 mm X 15 cm
FT-305.03	Ø 2.5 mm X 20 cm
FT-305.04	Ø 3.0 mm X 15 cm
FT-305.05	Ø 3.0 mm X 20 cm
FT-305.06	Ø 3.0 mm X 25 cm

Power Cannula

FT-306.01	Ø 2 mm X 15 cm
FT-306.02	Ø 2 mm X 20 cm
FT-306.03	Ø 2 mm X 25 cm
FT-306.04	Ø 2 mm X 30 cm
FT-306.05	Ø 2 mm X 35 cm
FT-306.06	Ø 3 mm X 15 cm
FT-306.07	Ø 3 mm X 20 cm
FT-306.08	Ø 3 mm X 25 cm
FT-306.09	Ø 3 mm X 30 cm
FT-306.10	Ø 3 mm X 35 cm
FT-306.11	Ø 4 mm X 15 cm
FT-306.12	Ø 4 mm X 20 cm
FT-306.13	Ø 4 mm X 25 cm
FT-306.14	Ø 4 mm X 30 cm
FT-306.15	Ø 4 mm X 35 cm
FT-306.16	Ø 5 mm X 15 cm
FT-306.17	Ø 5 mm X 20 cm
FT-306.18	Ø 5 mm X 25 cm
FT-306.19	Ø 5 mm X 30 cm
FT-306.20	Ø 5 mm X 35 cm
FT-306.21	Ø 6 mm X 15 cm
FT-306.22	Ø 6 mm X 20 cm
FT-306.23	Ø 6 mm X 25 cm
FT-306.24	Ø 6 mm X 30 cm
FT-306.25	Ø 6 mm X 35 cm



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula



Multiport Speed Harvester- One End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

60ccTommehub

FT-301.11	Ø 2.5 mm X 15 cm
FT-301.12	Ø 2.5 mm X 20 cm
FT-301.13	Ø 3.0 mm X 15 cm
FT-301.14	Ø 3.0 mm X 20 cm
FT-301.15	Ø 3.0 mm X 25 cm
FT-301.16	Ø 4.0 mm X 25 cm
FT-301.17	Ø 4.0 mm X 30 cm
FT-301.18	Ø 4.0 mm X 35 cm
FT-301.19	Ø 4.5 mm X 30 cm
FT-301.20	Ø 4.5 mm X 35 cm

Treaded Hub

FT-302.11	Ø 3.0 mm X 15 cm
FT-302.12	Ø 3.0 mm X 25 cm
FT-302.13	Ø 4.0 mm X 15 cm
FT-302.14	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-303.11	Ø 2.0 mm X 15 cm
FT-303.12	Ø 2.5 mm X 15 cm
FT-303.13	Ø 2.5 mm X 20 cm
FT-303.14	Ø 3.0 mm X 15 cm
FT-303.15	Ø 3.0 mm X 20 cm
FT-303.16	Ø 3.0 mm X 25 cm

Luer lock

FT-304.11	Ø 2.0 mm X 15 cm
FT-304.12	Ø 2.5 mm X 15 cm
FT-304.13	Ø 2.5 mm X 20 cm
FT-304.14	Ø 3.0 mm X 15 cm
FT-304.15	Ø 3.0 mm X 20 cm
FT-304.16	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-305.11	Ø 2.0 mm X 15 cm
FT-305.12	Ø 2.5 mm X 15 cm
FT-305.13	Ø 2.5 mm X 20 cm
FT-305.14	Ø 3.0 mm X 15 cm
FT-305.15	Ø 3.0 mm X 20 cm
FT-305.16	Ø 3.0 mm X 25 cm

Power Cannula

FT-306.11	Ø 2 mm X 15 cm
FT-306.12	Ø 2 mm X 20 cm
FT-306.13	Ø 2 mm X 25 cm
FT-306.14	Ø 2 mm X 30 cm
FT-306.15	Ø 2 mm X 35 cm
FT-306.16	Ø 3 mm X 15 cm
FT-306.17	Ø 3 mm X 20 cm
FT-306.18	Ø 3 mm X 25 cm
FT-306.19	Ø 3 mm X 30 cm
FT-306.20	Ø 3 mm X 35 cm
FT-306.21	Ø 4 mm X 15 cm
FT-306.22	Ø 4 mm X 20 cm
FT-306.23	Ø 4 mm X 25 cm
FT-306.24	Ø 4 mm X 30 cm
FT-306.25	Ø 4 mm X 35 cm
FT-306.26	Ø 5 mm X 15 cm
FT-306.27	Ø 5 mm X 20 cm
FT-306.28	Ø 5 mm X 25 cm
FT-306.29	Ø 5 mm X 30 cm
FT-306.30	Ø 5 mm X 35 cm
FT-306.31	Ø 6 mm X 15 cm
FT-306.32	Ø 6 mm X 20 cm
FT-306.33	Ø 6 mm X 25 cm
FT-306.34	Ø 6 mm X 30 cm
FT-306.35	Ø 6 mm X 35 cm



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula



Multi Port Harvester - Both End Bevel

A sophisticated, precision-designed harvester that is truly Patient Friendly. The ports in the Sorensen Harvester are machined with alternating 60 degree cutting edges. This unique design is less traumatic than other multi-port harvesters and delivers a fine-grained, high fat content/purity yield that can easily pass through 1.2mm, 0.9mm and 0.7mm blunt injectors.

60ccTommeyhub

FT-301.21	Ø 2.5 mm X 15 cm
FT-301.22	Ø 2.5 mm X 20 cm
FT-301.23	Ø 3.0 mm X 15 cm
FT-301.24	Ø 3.0 mm X 20 cm
FT-301.25	Ø 3.0 mm X 25 cm
FT-301.26	Ø 4.0 mm X 25 cm
FT-301.27	Ø 4.0 mm X 30 cm
FT-301.28	Ø 4.0 mm X 35 cm
FT-301.29	Ø 4.5 mm X 30 cm
FT-301.30	Ø 4.5 mm X 35 cm

Treaded Hub

FT-302.21	Ø 3.0 mm X 15 cm
FT-302.22	Ø 3.0 mm X 25 cm
FT-302.23	Ø 4.0 mm X 15 cm
FT-302.24	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-303.21	Ø 2.0 mm X 15 cm
FT-303.22	Ø 2.5 mm X 15 cm
FT-303.23	Ø 2.5 mm X 20 cm
FT-303.24	Ø 3.0 mm X 15 cm
FT-303.25	Ø 3.0 mm X 20 cm
FT-303.26	Ø 3.0 mm X 25 cm

Luer lock

FT-304.21	Ø 2.0 mm X 15 cm
FT-304.22	Ø 2.5 mm X 15 cm
FT-304.23	Ø 2.5 mm X 20 cm
FT-304.24	Ø 3.0 mm X 15 cm
FT-304.25	Ø 3.0 mm X 20 cm
FT-304.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-305.21	Ø 2.0 mm X 15 cm
FT-305.22	Ø 2.5 mm X 15 cm
FT-305.23	Ø 2.5 mm X 20 cm
FT-305.24	Ø 3.0 mm X 15 cm
FT-305.25	Ø 3.0 mm X 20 cm
FT-305.26	Ø 3.0 mm X 25 cm

Power Cannula

FT-306.21	Ø 2 mm X 15 cm
FT-306.22	Ø 2 mm X 20 cm
FT-306.23	Ø 2 mm X 25 cm
FT-306.24	Ø 2 mm X 30 cm
FT-306.25	Ø 2 mm X 35 cm
FT-306.26	Ø 3 mm X 15 cm
FT-306.27	Ø 3 mm X 20 cm
FT-306.28	Ø 3 mm X 25 cm
FT-306.29	Ø 3 mm X 30 cm
FT-306.30	Ø 3 mm X 35 cm
FT-306.31	Ø 4 mm X 15 cm
FT-306.32	Ø 4 mm X 20 cm
FT-306.33	Ø 4 mm X 25 cm
FT-306.34	Ø 4 mm X 30 cm
FT-306.35	Ø 4 mm X 35 cm
FT-306.36	Ø 5 mm X 15 cm
FT-306.37	Ø 5 mm X 20 cm
FT-306.38	Ø 5 mm X 25 cm
FT-306.39	Ø 5 mm X 30 cm
FT-306.40	Ø 5 mm X 35 cm
FT-306.41	Ø 6 mm X 15 cm
FT-306.42	Ø 6 mm X 20 cm
FT-306.43	Ø 6 mm X 25 cm
FT-306.44	Ø 6 mm X 30 cm
FT-306.45	Ø 6 mm X 35 cm



Injector Cannulas

(Micro Injectors)



Luer lock Luer lock (Stainless Steel) Ultra Luer lock Treaded 60cc Tommey hub Power Cannula

Sharp Extractor/Injector



This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

Ultra Luer lock

FT-403.01	Ø 1.0 mm X 5 cm
FT-403.02	Ø 1.0 mm X 7 cm
FT-403.03	Ø 2.0 mm X 15 cm
FT-403.04	Ø 2.5 mm X 15 cm
FT-403.05	Ø 3.0 mm X 15 cm

Luer lock

FT-404.01	Ø 1.0 mm X 5 cm
FT-404.02	Ø 1.0 mm X 7 cm
FT-404.03	Ø 2.0 mm X 15 cm
FT-404.04	Ø 2.5 mm X 15 cm
FT-404.05	Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.01	FT-405.11	FT-405.21	FT-405.31	FT-405.41
1.4Ø	FT-405.02	FT-405.12	FT-405.22	FT-405.32	FT-405.42
1.8Ø	FT-405.03	FT-405.13	FT-405.23	FT-405.33	FT-405.43
2.0Ø	FT-405.04	FT-405.14	FT-405.24	FT-405.34	FT-405.44

Blunt Extractor/Injector



This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

Ultra Luer lock

FT-403.11	Ø 2.0 mm X 15 cm
FT-403.12	Ø 2.5 mm X 15 cm
FT-403.13	Ø 3.0 mm X 15 cm

Luer lock

FT-404.11	Ø 2.0 mm X 15 cm
FT-404.12	Ø 2.5 mm X 15 cm
FT-404.13	Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.51	FT-405.61	FT-405.71	FT-405.81	FT-405.91
1.4Ø	FT-405.52	FT-405.62	FT-405.72	FT-405.82	FT-405.92
1.8Ø	FT-405.53	FT-405.63	FT-405.73	FT-405.83	FT-405.93
2.0Ø	FT-405.54	FT-405.64	FT-405.74	FT-405.84	FT-405.94

Blunt Tip Injector



Rounded blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp instruments.

Ultra Luer lock


FT-403.21	Ø 1.0 mm X 4 cm
FT-403.22	Ø 1.0 mm X 5 cm
FT-403.23	Ø 1.0 mm X 7 cm
FT-403.24	Curved 12 cm
FT-403.25	Curved 15 cm
FT-403.26	Curved 20 cm

Luer lock

FT-404.21	Ø 1.0 mm X 4 cm
FT-404.22	Ø 1.0 mm X 5 cm
FT-404.23	Ø 1.0 mm X 7 cm
FT-404.24	Curved 12 cm
FT-404.25	Curved 15 cm
FT-404.26	Curved 20 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.101	FT-405.111	FT-405.121	FT-405.131	FT-405.141
1.4Ø	FT-405.102	FT-405.112	FT-405.122	FT-405.132	FT-405.142
1.8Ø	FT-405.103	FT-405.113	FT-405.123	FT-405.133	FT-405.143
2.0Ø	FT-405.104	FT-405.114	FT-405.124	FT-405.134	FT-405.144



 **Fig. 2B**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.201	FT-405.211	FT-405.221	FT-405.231	FT-405.241
1.4Ø	FT-405.202	FT-405.212	FT-405.222	FT-405.232	FT-405.242
1.8Ø	FT-405.203	FT-405.213	FT-405.223	FT-405.233	FT-405.243
2.0Ø	FT-405.204	FT-405.214	FT-405.224	FT-405.234	FT-405.244

 **Fig. 1C**


Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.251	FT-405.261	FT-405.271	FT-405.281	FT-405.291
1.4Ø	FT-405.252	FT-405.262	FT-405.272	FT-405.282	FT-405.292
1.8Ø	FT-405.253	FT-405.263	FT-405.273	FT-405.283	FT-405.293
2.0Ø	FT-405.254	FT-405.264	FT-405.274	FT-405.284	FT-405.294

 **Fig. 1P**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.301	FT-405.311	FT-405.321	FT-405.331	FT-405.341
1.4Ø	FT-405.302	FT-405.312	FT-405.322	FT-405.332	FT-405.342
1.8Ø	FT-405.303	FT-405.313	FT-405.323	FT-405.333	FT-405.343
2.0Ø	FT-405.304	FT-405.314	FT-405.324	FT-405.334	FT-405.344

 **Fig. 2E**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.351	FT-405.361	FT-405.371	FT-405.381	FT-405.391
1.4Ø	FT-405.352	FT-405.362	FT-405.372	FT-405.382	FT-405.392
1.8Ø	FT-405.353	FT-405.363	FT-405.373	FT-405.383	FT-405.393
2.0Ø	FT-405.354	FT-405.364	FT-405.374	FT-405.384	FT-405.394



Luer lock

Luer lock
(Stainless Steel)

Ultra Luer lock

Treaded

60cc Tommey hub

Power Cannula

Coleman Liposuction Cannula Fig. MIC I

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm	3.5cm
1.1Ø	FT-405.401	FT-405.411	FT-405.421	FT-405.431	FT-405.441	FT-405.451	FT-405.461	FT-405.471
1.4Ø	FT-405.402	FT-405.412	FT-405.422	FT-405.432	FT-405.442	FT-405.452	FT-405.462	FT-405.472
1.8Ø	FT-405.403	FT-405.413	FT-405.423	FT-405.433	FT-405.443	FT-405.453	FT-405.463	FT-405.473
2.0Ø	FT-405.404	FT-405.414	FT-405.424	FT-405.434	FT-405.444	FT-405.454	FT-405.464	FT-405.474
2.5Ø	FT-405.405	FT-405.415	FT-405.425	FT-405.435	FT-405.445	FT-405.455	FT-405.465	FT-405.475

Coleman Liposuction Cannula Fig. MIC II

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	FT-405.501	FT-405.511	FT-405.521	FT-405.531	FT-405.541	FT-405.551	FT-405.561
1.4Ø	FT-405.502	FT-405.512	FT-405.522	FT-405.532	FT-405.542	FT-405.552	FT-405.562
1.8Ø	FT-405.503	FT-405.513	FT-405.523	FT-405.533	FT-405.543	FT-405.553	FT-405.563
2.0Ø	FT-405.504	FT-405.514	FT-405.524	FT-405.534	FT-405.544	FT-405.554	FT-405.564
2.5Ø	FT-405.505	FT-405.515	FT-405.525	FT-405.535	FT-405.545	FT-405.555	FT-405.565



Coleman Liposuction Cannula Spoon Tip Fig. MIC III

This tip works well for needle extraction of soft tissue and as a soft tissue re injection needle.

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	FT-405.601	FT-405.611	FT-405.621	FT-405.631	FT-405.641	FT-405.651	FT-405.661
1.4Ø	FT-405.602	FT-405.612	FT-405.622	FT-405.632	FT-405.642	FT-405.652	FT-405.662
1.8Ø	FT-405.603	FT-405.613	FT-405.623	FT-405.633	FT-405.643	FT-405.653	FT-405.663
2.0Ø	FT-405.604	FT-405.614	FT-405.624	FT-405.634	FT-405.644	FT-405.654	FT-405.664
2.5Ø	FT-405.605	FT-405.615	FT-405.625	FT-405.635	FT-405.645	FT-405.655	FT-405.665



Little Injector

Open and more pointed blunt tip is very effective in dense and fibrous tissue.

Ultra Luer lock

FT-403.51	Ø 1.0 mm X 4 cm
FT-403.52	Ø 1.0 mm X 5 cm
FT-403.53	Ø 1.0 mm X 7 cm

Luer lock

FT-404.51	Ø 1.0 mm X 4 cm
FT-404.52	Ø 1.0 mm X 5 cm
FT-404.53	Ø 1.0 mm X 7 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	FT-405.701	FT-405.711	FT-405.721	FT-405.731	FT-405.741
1.4Ø	FT-405.702	FT-405.712	FT-405.722	FT-405.732	FT-405.742
1.8Ø	FT-405.703	FT-405.713	FT-405.723	FT-405.733	FT-405.743
2.0Ø	FT-405.704	FT-405.714	FT-405.724	FT-405.734	FT-405.744




Fig. 1B

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	FT-405.801	FT-405.811	FT-405.821	FT-405.831	FT-405.841	FT-405.851	FT-405.861
1.4Ø	FT-405.802	FT-405.812	FT-405.822	FT-405.832	FT-405.842	FT-405.852	FT-405.862
1.8Ø	FT-405.803	FT-405.813	FT-405.823	FT-405.833	FT-405.843	FT-405.853	FT-405.863
2.0Ø	FT-405.804	FT-405.814	FT-405.824	FT-405.834	FT-405.844	FT-405.854	FT-405.864
2.5Ø	FT-405.805	FT-405.815	FT-405.825	FT-405.835	FT-405.845	FT-405.855	FT-405.865

Curved Injectors



Curved Injectors

New Curved Injectors are designed especially for fat transfer procedures to the body. Their unique shape mimics the anatomical curvature of the female torso and, thus, easily produces the natural-looking results that patients desire. Curved Injection Cannulas ameliorate the fat transfer process by providing a direct and precise route to the injection site and by lessening the time required to correctly sculpt curves.



Speciality Cannulas

Frizkon International Speciality Cannula are useful tools for different to operation are such as face or Neck.



V-Dissector

This cannula is specifically designed for use in the superficial compartment of the truncal area. The cutting notch is very useful for undermining in facial procedures and in dense, fibrous tissue.

60cc Tommeyhub

FT-501.01	Ø 2.5 mm X 15 cm
FT-501.02	Ø 2.5 mm X 20 cm
FT-501.03	Ø 3.0 mm X 15 cm
FT-501.04	Ø 3.0 mm X 20 cm
FT-501.05	Ø 3.0 mm X 25 cm
FT-501.06	Ø 4.0 mm X 25 cm
FT-501.07	Ø 4.0 mm X 30 cm
FT-501.08	Ø 4.0 mm X 35 cm
FT-501.09	Ø 4.5 mm X 30 cm
FT-501.10	Ø 4.5 mm X 35 cm

Treaded Hub

FT-502.01	Ø 3.0 mm X 15 cm
FT-502.02	Ø 3.0 mm X 25 cm
FT-502.03	Ø 4.0 mm X 15 cm
FT-502.04	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-503.01	Ø 2.0 mm X 15 cm
FT-503.02	Ø 2.5 mm X 15 cm
FT-503.03	Ø 2.5 mm X 20 cm
FT-503.04	Ø 3.0 mm X 15 cm
FT-503.05	Ø 3.0 mm X 20 cm
FT-503.06	Ø 3.0 mm X 25 cm

Luer lock


FT-504.01	Ø 2.0 mm X 15 cm
FT-504.02	Ø 2.5 mm X 15 cm
FT-504.03	Ø 2.5 mm X 20 cm
FT-504.04	Ø 3.0 mm X 15 cm
FT-504.05	Ø 3.0 mm X 20 cm
FT-504.06	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-505.01	Ø 2.0 mm X 15 cm
FT-505.02	Ø 2.5 mm X 15 cm
FT-505.03	Ø 2.5 mm X 20 cm
FT-505.04	Ø 3.0 mm X 15 cm
FT-505.05	Ø 3.0 mm X 20 cm
FT-505.06	Ø 3.0 mm X 25 cm

Power Cannula

FT-506.01	Ø 2 mm X 15 cm
FT-506.02	Ø 2 mm X 20 cm
FT-506.03	Ø 2 mm X 25 cm
FT-506.04	Ø 2 mm X 30 cm
FT-506.05	Ø 2 mm X 35 cm
FT-506.06	Ø 3 mm X 15 cm
FT-506.07	Ø 3 mm X 20 cm
FT-506.08	Ø 3 mm X 25 cm
FT-506.09	Ø 3 mm X 30 cm
FT-506.10	Ø 3 mm X 35 cm
FT-506.11	Ø 4 mm X 15 cm
FT-506.12	Ø 4 mm X 20 cm
FT-506.13	Ø 4 mm X 25 cm
FT-506.14	Ø 4 mm X 30 cm
FT-506.15	Ø 4 mm X 35 cm
FT-506.16	Ø 5 mm X 15 cm
FT-506.17	Ø 5 mm X 20 cm
FT-506.18	Ø 5 mm X 25 cm
FT-506.19	Ø 5 mm X 30 cm
FT-506.20	Ø 5 mm X 35 cm
FT-506.21	Ø 6 mm X 15 cm
FT-506.22	Ø 6 mm X 20 cm
FT-506.23	Ø 6 mm X 25 cm
FT-506.24	Ø 6 mm X 30 cm
FT-506.25	Ø 6 mm X 35 cm



Flap

The Flap is primarily used as an undermining dissector, and in smaller diameters, for facial tissue removal. The Flap Cannula in a 2.1mm diameter is very useful for undermining the nasolabial fold prior to autologous tissue injections, as well as for tissue removal in this anatomic area.

60ccTommeyhub

FT-501.11	Ø 2.5 mm X 15 cm
FT-501.12	Ø 2.5 mm X 20 cm
FT-501.13	Ø 3.0 mm X 15 cm
FT-501.14	Ø 3.0 mm X 20 cm
FT-501.15	Ø 3.0 mm X 25 cm
FT-501.16	Ø 4.0 mm X 25 cm
FT-501.17	Ø 4.0 mm X 30 cm
FT-501.18	Ø 4.0 mm X 35 cm
FT-501.19	Ø 4.5 mm X 30 cm
FT-501.20	Ø 4.5 mm X 35 cm

Treaded Hub

FT-502.11	Ø 3.0 mm X 15 cm
FT-502.12	Ø 3.0 mm X 25 cm
FT-502.13	Ø 4.0 mm X 15 cm
FT-502.14	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-503.11	Ø 2.0 mm X 15 cm
FT-503.12	Ø 2.5 mm X 15 cm
FT-503.13	Ø 2.5 mm X 20 cm
FT-503.14	Ø 3.0 mm X 15 cm
FT-503.15	Ø 3.0 mm X 20 cm
FT-503.16	Ø 3.0 mm X 25 cm

Luer lock


FT-504.11	Ø 2.0 mm X 15 cm
FT-504.12	Ø 2.5 mm X 15 cm
FT-504.13	Ø 2.5 mm X 20 cm
FT-504.14	Ø 3.0 mm X 15 cm
FT-504.15	Ø 3.0 mm X 20 cm
FT-504.16	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-505.11	Ø 2.0 mm X 15 cm
FT-505.12	Ø 2.5 mm X 15 cm
FT-505.13	Ø 2.5 mm X 20 cm
FT-505.14	Ø 3.0 mm X 15 cm
FT-505.15	Ø 3.0 mm X 20 cm
FT-505.16	Ø 3.0 mm X 25 cm

Power Cannula

FT-506.51	Ø 2 mm X 15 cm
FT-506.52	Ø 2 mm X 20 cm
FT-506.53	Ø 2 mm X 25 cm
FT-506.54	Ø 2 mm X 30 cm
FT-506.55	Ø 2 mm X 35 cm
FT-506.56	Ø 3 mm X 15 cm
FT-506.57	Ø 3 mm X 20 cm
FT-506.58	Ø 3 mm X 25 cm
FT-506.59	Ø 3 mm X 30 cm
FT-506.60	Ø 3 mm X 35 cm
FT-506.61	Ø 4 mm X 15 cm
FT-506.62	Ø 4 mm X 20 cm
FT-506.63	Ø 4 mm X 25 cm
FT-506.64	Ø 4 mm X 30 cm
FT-506.65	Ø 4 mm X 35 cm
FT-506.66	Ø 5 mm X 15 cm
FT-506.67	Ø 5 mm X 20 cm
FT-506.68	Ø 5 mm X 25 cm
FT-506.69	Ø 5 mm X 30 cm
FT-506.70	Ø 5 mm X 35 cm
FT-506.71	Ø 6 mm X 15 cm
FT-506.72	Ø 6 mm X 20 cm
FT-506.73	Ø 6 mm X 25 cm
FT-506.74	Ø 6 mm X 30 cm
FT-506.75	Ø 6 mm X 35 cm



Field Flap

With its elevated leading edge, this cannula is designed to open and expose soft tissue for more precise suctioning along the mandible and sub-mental areas.

60cc Tommey hub

FT-501.21	Ø 2.5 mm X 15 cm
FT-501.22	Ø 2.5 mm X 20 cm
FT-501.23	Ø 3.0 mm X 15 cm
FT-501.24	Ø 3.0 mm X 20 cm
FT-501.25	Ø 3.0 mm X 25 cm
FT-501.26	Ø 4.0 mm X 25 cm
FT-501.27	Ø 4.0 mm X 30 cm
FT-501.28	Ø 4.0 mm X 35 cm
FT-501.29	Ø 4.5 mm X 30 cm
FT-501.30	Ø 4.5 mm X 35 cm

Treaded Hub

FT-502.21	Ø 3.0 mm X 15 cm
FT-502.22	Ø 3.0 mm X 25 cm
FT-502.23	Ø 4.0 mm X 15 cm
FT-502.24	Ø 4.0 mm X 25 cm

Ultra Luer lock

FT-503.21	Ø 2.0 mm X 15 cm
FT-503.22	Ø 2.5 mm X 15 cm
FT-503.23	Ø 2.5 mm X 20 cm
FT-503.24	Ø 3.0 mm X 15 cm
FT-503.25	Ø 3.0 mm X 20 cm
FT-503.26	Ø 3.0 mm X 25 cm

Luer lock

FT-504.21	Ø 2.0 mm X 15 cm
FT-504.22	Ø 2.5 mm X 15 cm
FT-504.23	Ø 2.5 mm X 20 cm
FT-504.24	Ø 3.0 mm X 15 cm
FT-504.25	Ø 3.0 mm X 20 cm
FT-504.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

FT-505.21	Ø 2.0 mm X 15 cm
FT-505.22	Ø 2.5 mm X 15 cm
FT-505.23	Ø 2.5 mm X 20 cm
FT-505.24	Ø 3.0 mm X 15 cm
FT-505.25	Ø 3.0 mm X 20 cm
FT-505.26	Ø 3.0 mm X 25 cm

Power Cannula

FT-506.101	Ø 2 mm X 15 cm
FT-506.102	Ø 2 mm X 20 cm
FT-506.103	Ø 2 mm X 25 cm
FT-506.104	Ø 2 mm X 30 cm
FT-506.105	Ø 2 mm X 35 cm
FT-506.106	Ø 3 mm X 15 cm
FT-506.107	Ø 3 mm X 20 cm
FT-506.108	Ø 3 mm X 25 cm
FT-506.109	Ø 3 mm X 30 cm
FT-506.110	Ø 3 mm X 35 cm
FT-506.111	Ø 4 mm X 15 cm
FT-506.112	Ø 4 mm X 20 cm
FT-506.113	Ø 4 mm X 25 cm
FT-506.114	Ø 4 mm X 30 cm
FT-506.115	Ø 4 mm X 35 cm
FT-506.116	Ø 5 mm X 15 cm
FT-506.117	Ø 5 mm X 20 cm
FT-506.118	Ø 5 mm X 25 cm
FT-506.119	Ø 5 mm X 30 cm
FT-506.120	Ø 5 mm X 35 cm
FT-506.121	Ø 6 mm X 15 cm
FT-506.122	Ø 6 mm X 20 cm
FT-506.123	Ø 6 mm X 25 cm
FT-506.124	Ø 6 mm X 30 cm
FT-506.125	Ø 6 mm X 35 cm



Ultra Luer Lock

The Super Luer Lock is on our Single-use instruments.



Tumescent Riltrator

This cannula is used specifically for the iRlusion of tumescent solution and for hydraulic dissection. The configuration of the ports allows even distribution of the fluid.

SL-601.01 Ø 2.0 mm X 12cm
SL-601.02 Ø 2.0 mm X 20cm



Tonnard Harvester

A multiport harvester with 1mm elevated sharp microports. The 20 port configuration allows for harvesting both on the forward and backward strokes.

SL-602.01 Ø 3.0 mm X 20cm
SL-602.02 Ø 2.4 mm X 20cm



Tonnard Harvester

Excellent for harvesting small aliquots of tissues.

SL-603.01 Ø 2.5 mm X 15cm
SL-603.02 Ø 2.0 mm X 15cm



Carraway Harvester

Spiral three port design effective for multi-plane tissue harvesting.

SL-604.01 Ø 2.5 mm X 15cm
SL-604.02 Ø 2.0 mm X 10cm
SL-604.03 Ø 2.0 mm X 15cm
SL-604.04 Ø 1.5 mm X 8cm
SL-604.05 Ø 1.5 mm X 10cm



Injector

Round blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp injectors.

SL-605.01 Ø 1.5 mm X 9cm
SL-605.02 Ø 1.5 mm X 7cm
SL-605.03 Ø 1.0 mm X 5cm
SL-605.04 Ø 1.0 mm X 4cm



Spoon Tip Injector

Blunt injector designed with a more open port for rapid injection of larger volumes - a torso and buttock workhorse.

SL-606.01 Ø 3.0 mm X 20cm
SL-606.02 Ø 2.0 mm X 15cm



Miller Speed Harvester

SL-607.01 Ø 2.0 mm X 15cm

Liposuction Accessories

Handles

Connectors

- a) Luer to Luer
- b) Luer to Handle

Transfer Adapters

- a) Luer to Luer 1.2 mm
- Luer to Luer 1.4 mm
- Luer to Luer 2.4 mm
- b) Anaerobic Transfer
- c) Micro to Nano

Snap Lock

- a) Standard Snap Locks
60cc, 20cc, 10cc


b) frizkon International Spiral Snap Lock

Fat Transfer Gun

Decanting Stands


Sterilization Boxed

HANDLES FOR CANNULAS




FT-701.01

Luer Lock Handle



FT-702.01
20cc

Adaptor Handle Style I



FT-703.01
20cc

Adaptor Handle Style II

HANDLES FOR CANNULAS

FT-704.01

with connection thread for cannulas hose connection



FT-705.01

with connection thread for cannulas hose connection



FT-706.01

with connection thread for cannulas hose connection



FT-707.01

with connection thread for cannulas hose connection



FT-750.01

Pocar

The Pocar is a percutaneous punch that is used to make a circumferential incision prior to infiltration and tissue extraction. It is designed to maintain a tight seal while minimizing scarring.



NanoTransfer

FT-751.01

NanoTransfer




FT-752.01

Mesh for NanoTransfer



Anaerobic Transfers

- FT-760.01** Luer 1.2mm
- FT-760.02** Luer 1.4mm
- FT-760.03** Luer 2.4mm



FT-771.01

Luer to Luer



Decant Stands

FT-772.02 60cc to Luer
 FT-773.03 Luer to Luer



SnapLoks

These syringe snap locks fit into the lunger flanges of the syringe. When the plunger is pulled out the lock snaps open on the barrel lip thus holding the syringe under vacuum.

FT-901.01 60cc
 FT-901.02 50cc
 FT-901.03 35cc
 FT-901.04 20cc
 FT-901.05 10cc





Universal Power Injector

FT-801.01

The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or a 60cc Monoject syringe.



Fat Syringe

FT-950.01

The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or a 60cc Monoject syringe.




Sterilization Box

FT-1000.01



Syringe Decanting and Injection

FT-707.01	60 cc
FT-707.02	50 cc
FT-707.03	30 cc
FT-707.04	20 cc
FT-707.05	10 cc



FRIZKON

INTERNATIONAL

Pakistan Address:

Ireland Address:

USA Address:

Address:

No 665 Naaz Street K LW Sialkot 51310

Email:

info@frizkontools.com

Contact:

+923314494637

+923325334864

+0523306091

Website:

www.frizkontools.com

Address:

35 Belmont Middle Gardiner Street

Dublin D01C447

Email:

info@frizkontools.com

Contact:

+353831052805

Website:

www.frizkontools.com

Address:

13626 Near Point Tomball TX77377

Email:

info@frizkontools.com

Contact:

+1 (832) 686-3125

Website:

www.frizkontools.com